

JUL 11 2005

K 051752

510(k) Summary – Elecsys® Troponin T STAT Assay

Introduction According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence

Submitter name, address, contact Roche Diagnostics
9115 Hague Rd
Indianapolis IN 46250
(317) 521-3723

Contact person: Theresa M. Ambrose

Date prepared: June 28, 2005

Device Name Proprietary name: Elecsys® Troponin T STAT Assay

Common name: Troponin T Test

Classification name: Immunoassay method, troponin subunit

Device Description The Elecsys® Troponin T STAT assay is a two step sandwich immunoassay with streptavidin micro particles and electrochemiluminescence detection, for the measurement of human TnT in serum or plasma .

Intended use Immunoassay for the in vitro quantitative determination of troponin T in human serum and plasma. Elecsys Troponin T can be used as an aid in the differential diagnosis of acute coronary syndrome to identify necrosis, e.g., acute myocardial infarction. The test is further indicated for the risk stratification of patients presenting with acute coronary syndrome and for cardiac risk in patients with chronic renal failure. The test may also be useful for the selection of more intensive therapy and intervention in patients with elevated levels of cardiac Troponin T. The electrochemiluminescence immunoassay "ECLIA" is intended for use on the Roche Elecsys family of immunoassay analyzers.

Predicate Device We claim substantial equivalence to the Elecsys ® Troponin T STAT reagent currently marketed on the Elecsys 2010 and MODULAR ANALYTICS E170. (K040733).

Continued on next page

510(k) Summary – Elecsys® Troponin T STAT Assay, Continued

Substantial equivalency – device comparison The table below indicates the similarities and differences between the modified Elecsys® Elecsys Troponin T STAT and its predicate device (current Troponin T STAT, K040733).

Characteristic	Current Elecsys® Troponin T STAT (3 rd Gen.) (K040733)	Modified device Elecsys® Troponin T STAT (4 th Gen.)
Intended Use	<p>Immunoassay for the in vitro quantitative determination of troponin T in human serum and plasma. Elecsys Troponin T can be used as an aid in the differential diagnosis of acute coronary syndrome to identify necrosis, e.g., acute myocardial infarction. The test is further indicated for the risk stratification of patients presenting with acute coronary syndrome and for cardiac risk in patients with chronic renal failure. The test may also be useful for the selection of more intensive therapy and intervention in patients with elevated levels of cardiac Troponin T. The electrochemiluminescence immunoassay “ECLIA” is intended for use on the Roche Elecsys family of immunoassay analyzers.</p>	Same
Indications for Use	<p>Elecsys Troponin T can be used as an aid in the differential diagnosis of acute coronary syndrome to identify necrosis, e.g., acute myocardial infarction. The test is further indicated for the risk stratification of patients presenting with acute coronary syndrome and for cardiac risk in patients with chronic renal failure. The test may also be useful for the selection of more intensive therapy and intervention in patients with elevated levels of cardiac Troponin T.</p>	Same

Special 510(k): Device Modification – Elecsys® Troponin T STAT ,Continued

Assay principle	Electrochemiluminescent immunoassay	Same
Traceability/standardization	Standardized against the 2 nd generation Troponin T test	Standardized against the 3rd generation Elecsys Troponin T test; traceable to the 2 nd generation test
Calibration frequency	<p>Elecsys 2010</p> <ul style="list-style-type: none"> • After 1 month when using the same reagent lot • After 7 days when using the same reagent kit on the analyzer <p>Elecsys 1010</p> <ul style="list-style-type: none"> • With every reagent kit • After 7 days (20-25 °C) • After 3 days (25-32 °C) 	<p>Elecsys 2010</p> <ul style="list-style-type: none"> • After 1 month when using the same reagent lot • After 7 days when using the same reagent kit on the analyzer <p>Elecsys 1010</p> <ul style="list-style-type: none"> • With every reagent kit • After 7 days (20-25 °C)
Sample type	Human serum, K3-EDTA and Na-citrate plasma.	Human serum, K2- and K3-EDTA, Li-heparin, and Na-citrate plasma
Reagent stability	<p>Unopened</p> <ul style="list-style-type: none"> • Up to the stated expiration date at 2-8 °C <p>After opening</p> <ul style="list-style-type: none"> • 12 weeks at 2-8 °C • 8 weeks on Elecsys 2010 • 8 weeks on Elecsys 1010 (20-25 °C; up to 20 hours opened in total) 	Same
Calibrator	Elecsys Troponin T STAT CalSet	Same
Controls	Elecsys PreciControl Troponin T or Elecsys PreciControl Cardiac.	Elecsys PreciControl Troponin T
Duration of assay	9 minutes	Same

Special 510(k): Device Modification – Elecsys® Troponin T STAT ,Continued

Measuring range	0.010 – 25.00 ng/mL	Same
Precision	<p>Within-run (human serum)</p> <ul style="list-style-type: none"> • 1.1% CV at 0.47 ng/mL • 1.1% CV at 2.63 ng/mL • 1.4% CV at 11.5 ng/mL <p>Within-run (PeciControl)</p> <ul style="list-style-type: none"> • 4.2% CV at 0.10 ng/mL • 3.0% CV at 5.07 ng/mL <p>Total (human serum)</p> <ul style="list-style-type: none"> • 5.8% CV at 0.47 ng/mL • 5.4% CV at 2.63 ng/mL • 5.7% CV at 11.5 ng/mL <p>Total (PeciControl)</p> <ul style="list-style-type: none"> • 9.3% CV at 0.10 ng/mL • 6.0% CV at 5.07 ng/mL 	<p>Within-run (human serum)</p> <ul style="list-style-type: none"> • 4.5% CV at 0.047 ng/mL • 2.0% CV at 0.652 ng/mL • 2.9% CV at 6.08 ng/mL <p>Within-run (PeciControl)</p> <ul style="list-style-type: none"> • 2.2% CV at 0.137 ng/mL • 2.5% CV at 2.89 ng/mL <p>Total (human serum)</p> <ul style="list-style-type: none"> • 6.2% CV at 0.047 ng/mL • 4.6 % CV at 0.652 ng/mL • 5.6% CV at 6.08 ng/mL <p>Within-run (PeciControl)</p> <ul style="list-style-type: none"> • 3.5% CV at 0.137 ng/mL • 4.7% CV at 2.89 ng/mL
Concentration at 10% CV:	0.03 ng/mL	Same
Hook effect	No hook effect up to 400 ng/mL	Same
Analytical sensitivity	Lower detection limit: 0.01 ng/mL Lowest concentration giving 10 % CV: 0.03 ng/mL	Same

Special 510(k): Device Modification – Elecsys® Troponin T STAT, Continued

Limitations – interference	No interference from <ul style="list-style-type: none"> • icterus up to 27 mg/dL bilirubin • hemolysis up to 0.1 g/dL, • Lipemia up to 1500 mg/dL Intralipid • Biotin up to 50 ng/mL • Rheumatoid factor up to 2000 U/mL 	Same
	Falsely depressed results are obtained when using samples with higher hemoglobin concentrations.	Same
	Plasma samples collected with heparin or oxalate/fluoride revealed sample-dependent low TnT values compared to results obtained on serum samples.	Plasma samples collected using tubes containing oxalate/fluoride revealed sample-dependent low TnT values when compared to results obtained on serum samples.
	Contains additives to minimize the effects of interference due to <ul style="list-style-type: none"> • Monoclonal mouse antibodies • Antibodies to streptavidin 	Same
	Extremely high titers of antibodies to ruthenium can cause interference.	Same
	Results should be assessed in conjunction with the patient’s medical history, clinical examination, and other findings .	Same



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

JUL 11 2005

Theresa M. Ambrose, Ph.D., RAC
Regulatory Principal
Roche Diagnostics
9115 Hague Road
Indianapolis, IN 46250

Re: k051752
Trade/Device Name: Roche Elecsys® Troponin T STAT
Regulation Number: 21 CFR 862.1215
Regulation Name: Creatine phosphokinase/creatin kinase or isoenzymes test system
Regulatory Class: Class II
Product Code: MMI
Dated: June 28, 2005
Received: June 29, 2005

Dear Dr. Ambrose:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Carol C. Benson, M.A.
Acting Director
Division of Chemistry and Toxicology
Office of *In Vitro* Diagnostic Device
Evaluation and Safety
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K051752

Device Name: Roche Elecsys® Troponin T STAT

Indications For Use:

Immunoassay for the in vitro quantitative determination of troponin T in human serum and plasma. Elecsys Troponin T can be used as an aid in the differential diagnosis of acute coronary syndrome to identify necrosis, e.g., acute myocardial infarction. The test is further indicated for the risk stratification of patients presenting with acute coronary syndrome and for cardiac risk in patients with chronic renal failure. The test may also be useful for the selection of more intensive therapy and intervention in patients with elevated levels of cardiac Troponin T. The electrochemiluminescence immunoassay "ECLIA" is intended for use on the Roche Elecsys family of immunoassay analyzers.

Prescription Use XXX
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)


Division Sign-Off

Office of In Vitro Diagnostic Device Page 1 of 1
Evaluation and Safety

Confidential

510(k) K051752

22