

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)

Iso-Compass User Manual

Isotope Data Reduction Software

(Ver. 1.0)

Wen Zhang

State Key Laboratory of Geological Processes and Mineral Resources

China University of Geosciences (Wuhan)

No.388 Lumo Road, Wuhan city, Hubei province, P.R. China

E-mail: tuyaken@hotmail.com



Contents

1. Software introduction	3
1.1 Iso-Compass	3
1.2 Software features	3
1.3 Operating environment	3
2. Software installation	5
3. Software workflow and functions	8
3.1 Software workflow	8
3.2 Select initial file information	10
3.3 Create an analysis method template	12
3.3.1 Read data	12
3.3.2 "Formula" module	14
3.3.3 "Normalization" module	17
3.3.4 "Export" module	19
3.4 Data reduction	20
3.5 Saving data and exporting reports	22
4. Laser-MC-ICP-MS data reduction	25
4.1 Data preparation	25
4.2 Data import	25
4.3 Data reduction	26
4.4 Data export	27

1. Software introduction

1.1 Iso-Compass

Iso-Compass is an isotope composition data reduction software developed by the LA-ICP-MS Division of **the State Key Laboratory of Geological Processes and Mineral Resources, China University of Geosciences (Wuhan)**. The software is mainly used for processing raw isotope data in a batch and pattern manner, with the aim of processing the isotope data rapidly and flexibly.

1.2 Software features

1. Reads MC-ICP-MS raw data simply and conveniently; combines graphs and tables to display data intuitively and accurately.
2. Flexibly writes data correction formulas, including instrument mass fraction correction (such as internal standard correction, SSB correction, and external element correction) and interference correction.
3. Satisfies solution MC-ICP-MS data reduction and laser MC-ICP-MS data reduction to achieve batch processing of high-precision isotopic composition absolute and delta data.
4. Quickly formulates the output report format to achieve a formatted report.
5. The software is based on the .Net framework and it can be used independently on the Windows operating system, without installing other commercial software.

1.3 Operating environment

Hardware requirements

- Recommended minimum configuration: Pentium 1 GHz or above, 512 MB RAM or above
- Minimum disk space: x86 – 850 MB; x64 – 2 GB

Applicable software

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)

- Windows XP SP3
- Windows Vista SP1 or above
- Windows Server 2008 (not supported on Server Core roles)
- Windows Server 2008 R2 (not supported on Server Core roles)
- Windows 7 SP1 operating system and above
- Windows Server 2008 R2 SP1 operating system and above

2. Software installation

The IsoCompass installation package contains the following files. Click "setup.exe" to start the installation.

名称	修改日期	类型	大小
DotNetFX46	2019/1/16 22:38	文件夹	
setup.exe	2019/1/16 22:38	应用程序	772 KB
Setup.msi	2019/1/16 22:38	Windows Install...	1,699 KB

Figure 2.1 IsoCompass installation package

When the dialog box displayed in Fig. 2.2 pops up, click "Next."




Figure 2.2 Installation process

Click "Browse" and select the software installation path.




Figure 2.3 Installation process

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)

After selecting the installation path, click "Next." Once the software installation is complete, click "Close" to exit.




Figure 2.4 Software installation complete

This software comes with the Microsoft .NET Framework 4 installation package. If the installation system finds that the pre-installed Microsoft .NET Framework 4 is missing, the installation system will first install Microsoft .NET Framework 4, and then install the Iso-Compass software.

The Iso-Compass software following installation is illustrated in Fig. 2.5.




Figure 2.5 Iso-Compass folder

Once the software installation is complete, you need to register. Firstly,

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)

click on the "IsoCompass.exe" file. Please open the file in “**administrator mode.**” If you are using the software for the first time, the following dialog box will pop up:




Figure 2.6 Verify software license

Please click "Apply License." The software will generate a document, "Requestcode.txt." Please send this document to the e-mail address "tuyaken@hotmail.com" to apply for a software license from the author. Please search for "Requestcode.txt" in the software installation folder.

After obtaining a valid software license, you can directly find the path where the software license is located. After confirming the license path, the software will prompt the validity period of the license. Please click "OK" to complete the verification.





Figure 2.7 Software license verified

If the software license has expired, please send an e-mail to "tuyaken@hotmail.com" to obtain a new license.

3. Software workflow and functions

3.1 Software workflow



3.2 Select initial file information

After clicking the "IsoCompass.exe" file, the welcome screen will appear, followed by the initial dialog box of the software. In the initial dialog, you need to set "Workspace (data storage path)," "Folder (raw data storage path)," "File Type (raw data file format)," "SampleList (raw data list/optional)," and "Approach (Analysis method template)."




Figure 3.1 Iso-Compass data import interface

Workspace (data save path)

This is used to store all of the data output by the Iso-Compass software, including data reports saved in the output report and data reduction.

Folder (raw data storage path)

The folder in which the data to be processed are located.




Figure 3.2 Iso-Compass data import interface


Iso-Compass Isotope Data Reduction Software (Ver. 1.0)

File Type (raw data file format)

In the folder where the data to be processed exist, if there are different document types, you can select the file format of the data to be processed. For example, the raw data files exported by Neptune Plus instruments are in the ".exp" format.

SampleList (raw data list/optional)

As Iso-Compass can read all data in the specified file format in "Folder (raw data storage path)," Iso-Compass can read data without relying on "SampleList." However, the data only have the analysis test name, and no sample name. If you need the sample name for external standard calibration, you need to prepare a sample list (Excel format). The sample list needs to be recorded in the specified format, as illustrated in Fig. 3.3. Iso-Compass will depend on the file name in "Folder (raw data storage path)" to find the corresponding sample name in "SampleList."



	A	B	C	D	E	F
1	YANG	JAN26B001	BHVO-2G			
2	YANG	JAN26B002	BHVO-2G			
3	YANG	JAN26B003	BCR-2G			
4	YANG	JAN26B004	HNB-08			
5	YANG	JAN26B005	HNB-08			
6	YANG	JAN26B006	T16-45-01			
7	YANG	JAN26B007	T16-45-02			
8	YANG	JAN26B008	T16-45-03			
9	YANG	JAN26B009	T16-45-04			
10	YANG	JAN26B010	T16-45-05			
11	YANG	JAN26B011	HNB-08			
12	YANG	JAN26B012	HNB-08			
13	YANG	JAN26B013	BHVO-2G			
14	YANG	JAN26B014	BHVO-2G			
15	YANG	JAN26B015	T16-45-06			

Figure 3.3 List file writing format

"Set as default"

This sets the parameters set above to the default values.

After setting "Workspace (data storage path)," "Folder (raw data storage path)," and "File Type (raw data file format)," you can select an edited analysis method template in "Approach (analysis method template)" and then click "OK" to start data reduction. You can also click "Preview/New" in "Approach" to view or modify the edited analysis method template.

If no suitable analysis method template is available, you can click "Preview/New" to create a new analysis method template (see following section).

You can click "Load Session" to read the previously saved session data.

3.3 Create analysis method template

3.3.1 Read data

Click "Preview/New" in "Approach" and an initial analysis method template will appear.




Figure 3.1 Initial analysis method template interface

First click "New" and enter the name of the new analysis method.

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)




Figure 3.2 Enter name of new analysis method

Thereafter, you need to set the "row" and "column" path of the file data to be processed. You need to open a raw file of the data to be processed, as with the Sr isotope raw data file displayed below:


R152														
A	B	C	D	E	F	G	H	I	J	K	L	M		
Triton Analysis Data Report														
Run-Number-View														
Filename: F:\2019_LASER\MC19_025_YANGXIAOSONG\JAN26B01.dat														
Date: 2019/1/28														
Sample ID:														
Method Name:														
Wheel ID:														
Run number: 0														
Analysis date: 2019/1/28														
Analysis time: 10:59:23														
Operator:														
Instrument:														
Comment:														
Cycle	Time	83Kr	167Er++	84Sr	85Rb	85.9	173Yb++	87Sr	88Sr					
16	1 14:39:55.4	1.83E-04	4.57E-05	1.06E-03	1.07E-05	3.10E-04	-7.01E-05	1.33E-04	-8.97E-05					
17	2 14:39:55.9	2.13E-04	-5.38E-05	1.06E-03	-5.66E-06	2.58E-04	-5.56E-05	3.15E-05	5.89E-05					
18	3 14:39:56.5	2.63E-04	2.25E-05	1.06E-03	-4.43E-05	3.41E-04	-5.87E-05	-5.31E-05	5.77E-05					
19	4 14:39:57.0	1.85E-04	1.78E-05	9.78E-04	-5.89E-06	3.46E-04	-6.78E-05	-4.84E-05	1.88E-05					
20	5 14:39:57.5	1.95E-04	-3.55E-06	1.06E-03	-6.72E-05	2.71E-04	-1.60E-05	3.34E-05	1.36E-05					
21	6 14:39:58.0	1.62E-04	2.58E-05	1.12E-03	-4.52E-05	4.80E-04	-2.09E-05	7.61E-07	-5.90E-06					
22	7 14:39:58.6	2.42E-04	8.74E-06	9.92E-04	1.58E-04	2.81E-04	-5.45E-05	-2.12E-05	3.09E-05					
23	8 14:39:59.1	2.31E-04	2.39E-05	1.03E-03	1.14E-04	3.57E-04	-6.72E-06	-6.18E-05	3.35E-05					
24	9 14:39:59.6	3.01E-04	-2.82E-05	1.02E-03	1.06E-04	3.93E-04	-3.66E-05	-2.78E-05	3.75E-05					
25	10 14:40:00.1	2.57E-04	-2.38E-05	1.05E-03	4.41E-06	3.26E-04	-4.73E-05	-2.29E-05	4.06E-05					
26	11 14:40:00.7	1.88E-04	-2.30E-05	1.14E-03	7.58E-05	3.96E-04	-4.62E-05	-5.29E-05	9.83E-07					
27	12 14:40:01.2	1.91E-04	6.69E-06	1.07E-03	-1.08E-05	3.98E-04	-4.83E-05	-1.12E-04	-1.63E-06					
28	13 14:40:01.7	2.54E-04	1.17E-05	1.01E-03	-1.27E-05	2.81E-04	-9.05E-05	-4.21E-07	-8.12E-05					
29	14 14:40:02.3	2.10E-04	2.67E-05	9.91E-04	1.05E-05	3.29E-04	-6.19E-05	9.48E-05	1.08E-04					
30	15 14:40:02.8	2.18E-04	-5.29E-05	9.85E-04	-1.43E-05	3.32E-04	-4.55E-05	-3.99E-05	1.78E-05					
31	16 14:40:03.3	2.35E-04	1.01E-05	1.07E-03	2.36E-05	2.74E-04	-4.06E-05	-9.16E-06	-6.38E-06					
32	17 14:40:03.8	1.83E-04	7.18E-06	1.01E-03	3.13E-05	3.38E-04	-1.83E-05	3.45E-05	9.06E-06					
33	18 14:40:04.4	1.81E-04	-1.84E-05	1.03E-03	7.53E-05	2.52E-04	-2.83E-05	-4.58E-05	-4.70E-05					
34	19 14:40:04.9	2.00E-04	-3.57E-05	1.12E-03	5.92E-05	2.98E-04	-9.55E-06	-3.75E-05	1.43E-05					
35	20 14:40:05.4	1.92E-04	2.55E-05	1.04E-03	9.32E-06	3.27E-04	-6.24E-05	-7.11E-05	7.20E-05					
36	21 14:40:05.9	2.02E-04	1.75E-05	9.43E-04	6.43E-05	3.28E-04	-4.90E-05	-2.48E-05	-5.62E-05					
37	22 14:40:06.5	2.57E-04	7.44E-06	1.03E-03	7.92E-06	3.07E-04	-5.28E-05	-7.22E-05	1.19E-05					
38	23 14:40:07.0	2.61E-04	3.43E-05	9.87E-04	3.83E-05	3.96E-04	-2.86E-05	-1.13E-05	7.10E-05					
39	24 14:40:07.5	1.94E-04	4.27E-05	1.06E-03	-1.97E-05	3.47E-04	-3.31E-05	1.94E-06	-1.04E-05					
40	25 14:40:08.0	2.04E-04	-3.17E-05	9.67E-04	-2.79E-05	4.05E-04	-2.18E-05	-2.71E-05	3.07E-05					
41	26 14:40:08.6	2.02E-04	-6.55E-05	1.03E-03	-3.14E-05	2.45E-04	-4.87E-05	6.38E-05	4.42E-05					
42	27 14:40:09.1	2.29E-04	-3.40E-05	9.76E-04	-2.48E-05	3.50E-04	-7.10E-05	-1.60E-06	2.02E-05					
43	28 14:40:09.6	2.00E-04	1.10E-05	1.09E-03	3.74E-05	2.39E-04	-6.84E-05	-1.28E-04	9.53E-06					
44	29 14:40:10.2	2.50E-04	-9.65E-06	1.02E-03	-1.81E-05	4.43E-04	-2.21E-05	5.49E-06	-5.62E-05					
45	30 14:40:10.7	2.19E-04	-5.22E-05	1.06E-03	4.09E-05	3.44E-04	-1.85E-05	-3.33E-05	3.40E-05					
46	31 14:40:11.2	2.26E-04	-5.29E-05	1.10E-03	-3.32E-06	2.60E-04	-8.96E-05	7.30E-05	5.97E-06					
47	32 14:40:11.7	2.41E-04	1.59E-05	8.98E-04	-1.76E-05	1.59E-04	-3.76E-05	1.35E-05	5.02E-06					
48	33 14:40:12.3	1.68E-04	3.05E-05	9.35E-04	6.60E-05	3.24E-04	-3.98E-06	4.49E-05	-7.36E-05					
49	34 14:40:12.8	1.93E-04	2.32E-05	1.04E-03	1.82E-05	3.46E-04	-6.19E-05	1.38E-05	-4.98E-05					
50	35 14:40:13.3	1.54E-04	1.84E-05	1.09E-03	4.44E-05	3.91E-04	-3.27E-05	-1.20E-04	-3.27E-05					
51	36 14:40:13.8	2.15E-04	-8.23E-06	9.75E-04	7.32E-05	3.54E-04	-4.51E-05	4.66E-05	-1.68E-05					
52	37 14:40:14.4	2.02E-04	-1.50E-05	9.65E-04	-7.42E-07	2.74E-04	5.42E-06	-5.38E-06	-4.44E-05					

Figure 3.3 Sr isotope raw data file of Neptune Plus (in ".exp" format)

You should find the start and end numbers of the row and column of the data to read them. For example, the data start at row 16 and end at row 175, and start at column C (column 3) and end at column J (column 10).

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)


Enter the data into "Row" and "Column."



The screenshot shows a software window with a title bar. Below the title bar, there are two input fields. The first is labeled "Row" and contains the number "16". To its right is a minus sign and another input field containing "175". The second is labeled "Column" and contains the number "3". To its right is a minus sign and another input field containing "10".

Figure 3.4 Enter row and column numbers of data to read them

You can design the background and signal interval in the imported data in advance. For example, the 1st to 45th data points are background intervals, while the 65th to 140th data points are signal intervals.



The screenshot shows a software window with a title bar. Below the title bar, there are four buttons: "Save", "Save As...", "New", and "Delete". Below these buttons, there are two input fields. The first is labeled "Background" and contains the number "1". To its right is a minus sign and another input field containing "45". The second is labeled "Signal" and contains the number "65". To its right is a minus sign and another input field containing "140".

Figure 3.5 Set signal acquisition interval of data

3.3.2 "Formula" module

After setting the "row" and "column" locations of the data, "Raw Data" will display the name of the isotope to be tested for each data column. Note: "r_1, r_2, r_x" are used as the code of the corresponding isotope, which is used for subsequent formula editing.

Check the small box under "Bkg Substraction" and the ticked "Bkg Substraction" box to which it refers will perform background correction. The background deduction principle of the Iso-Compass software is as follows: firstly, select the interval of the background signal through "Background," then remove the 2SD outlier data from the background interval and calculate the average of the remaining data as the "background value." Data with the ticked "Bkg Substraction" box will subtract "background value" from the original data.

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)

Formula			Normalization	Export
Raw Data				
Data Name	APR29B001.static	Bkg Subtraction		
r_1	83Kr	<input checked="" type="checkbox"/>		
r_2	167Er++	<input checked="" type="checkbox"/>		
r_3	84Sr	<input checked="" type="checkbox"/>		
r_4	85Rb	<input checked="" type="checkbox"/>		
r_5	85.90	<input checked="" type="checkbox"/>		
r_6	173Yb++	<input checked="" type="checkbox"/>		
r_7	87Sr	<input checked="" type="checkbox"/>		
r_8	88Sr	<input checked="" type="checkbox"/>		

Figure 3.6 Original signal reading and background value subtraction


You can edit formulas in "Formula," including calculation methods such as addition, subtraction, multiplication, division, and exponent. As indicated in the figure below, "f1" refers to a formula code, which can be used for formula calculation. "Name" refers to the actual formula name and it cannot be used in calculations, but it is indicated as the formula name in the following diagrams. You can enter specific formulas in "Expression."

Formula					
Name	Expression	Mean	2SE	Count	
f1	84Sr	$r_3 - r_1 * c1 * (57/11.49) - r_2 * (26.978/22.869)$	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
f2	85Rb	$r_4 - r_2 * (14.91/22.869) - r_6 * (2.982/16.103)$	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
f3	87Sr-Rb	$r_7 - r_6 * (32.26/16.103)$	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
f4	86Sr	$r_5 - r_1 * c1 * (17.3/11.49) - r_6 * (21.68/16.103)$	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
f5	88Sr	$r_8 - r_6 * (12.996/16.103)$	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
f6	F-Sr	$\ln(8.37521 / ((f5/f4) / \ln(87.90562/85.90927)))$	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
f7	F-Rb	$\ln(0.38571 / (((f3 - f4 * c2) / ((87/86)^{f6})) / f2)) / \ln(87/85)$	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
f8	F-Rb/F-Sr	$f7/f6$	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
f9	87Sr	$f3 - f2 * 0.38571 / (((87/85)^{f6 * c3}))$	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
f11	87Sr/86Sr	$f9/f4 * ((87/86)^{f6})$	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
f1	84Sr/86Sr	$f1/f4 * ((84/86)^{f6})$	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
f1	83Kr	r_1	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
f1	Rb/Sr	$(f2/f5) * 1.144$	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Figure 3.7 Data reduction and editing formulas

For example, the original signal of ⁸⁴Sr is (r_3), and the interference from ⁸³Kr+(r_1) and ⁶⁷Er++ (r_2) needs to be subtracted. Then, the correction formula is " $r_3 - r_1 * c1 * (57/11.49) - r_2 * (26.978/22.869)$," where "c1" is a constant parameter set in "Constants." "Constants" can be used to record certain commonly used constant parameters, but the value must be a number and not a formula.

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)



	Name	Value	
	c1	F-Kr	1
	c2	87Sr/86Sr	0.705003
	c3	f-Rb/f-Sr	1.085

Figure 3.8 Designing constant terms in formula editing

Note: The data in "Formula" are calculated based on the data in each "cycle" (the smallest data unit), following which the average value and relative standard error (2SE) are provided.

Calculations supported by "Formula":

Addition: "+"

Subtraction: "-"


Multiplication: "*"

Division: "/"

Index: "^"

Natural logarithm: "ln()"

Iso-Compass provides a special constant calculation mode, "Average Parameter." The calculation module supports entering a calculation formula to obtain a constant. That is, the user can enter a calculation formula in this module; for example, the calculation formula for calculating the Hf isotope fractionation factor in Hf isotope analysis, and then the formula outputs the average value of the fractionation factor. The user can opt to use the average value of the fractionation factor to correct the mass fractionation or use the fractionation factor in "Formula" to perform a one-to-one mass fractionation correction.



	Name	Value	Reference	
	a1		No External	

Figure 3.9 Designing variable constants in formula editing

The "Average Parameter" function can also be used. You need to use a calibration parameter in the reference material to correct the data in the actual analysis sample. For example, if you need to correct the interference of Rb on the Sr isotope, you can analyze a reference substance containing Rb before analyzing the actual sample, calculate the Rb isotope fractionation factor from

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)

the reference substance, and then use this isotope fractionation factor for the actual sample. You simply need to select the desired "Reference" name in "Average Parameter" (such as NIST610). This indicates that the constant "a1" is calculated by all selected NIST610 in this data reduction, and this constant can subsequently be used in the "Formula" formula.




Figure 3.10 Designing variable constants in formula editing

3.3.3 "Normalization" module

The "Normalization" module is used to implement external standard correction. External standard calibration first requires the user to provide an accurate sample name. A reference material data file "ReferenceData.exp" (C:\IsoCompass\config) is provided in Iso-Compass. The user can add the required reference substance name to the database, and enter the isotope name and accurate isotope ratio.

	A	B	C	D	E	F	G	H
1		208Pb/206Pb	207Pb/206Pb	208Pb/204Pb	207Pb/204Pb	206Pb/204Pb	87Sr/86Sr	84Sr/86Sr
2	NIST610	2.169	0.91	36.991	15.515	17.052		
3	NIST612	2.165	0.907	37.02	15.516	17.099		
4	BCR-2G	2.063	0.832	38.709	15.614	18.761	0.705003	0.05657
5	BHVO-2G			38.273	15.591	18.778	0.703469	0.05657
6	BIR-2G						0.703105	0.05657
7	NBS981	2.168	0.915	36.726	15.5	16.942		
8	NBS987						0.710241	0.05657

Figure 3.11 Reference material data file "ReferenceData.exp"

You can find the names of the reference materials in the database in "Reference Materials" of "Normalization" in the software and then select the required reference materials.

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)




Figure 3.12 Selecting externally calibrated reference substance

In "Corrected Isotopic Ratio," select the isotope ratio of the reference material in the "Reference Isotop" column and then select the name of the formula written in the "Formula" module in the actual analysis in "Measured Data". As indicated in the figure below, $^{87}\text{Sr}/^{86}\text{Sr}$ is provided to correct the reference value of the reference substance.




Figure 3.13 Selecting isotope ratio for external calibration

After selecting the ratio of the reference material for the isotope to be corrected, select the correction scheme at the bottom:

No Correction: No external standard correction.

Mean Correction: Use the reference material before and after to calculate

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)

the correction factor, and take the average of the correction factor to correct the unknown sample.

Interpolation Correction: Calculate the average value of the front and back reference material separately. It is assumed that the isotope fractionation changes between the front and back reference materials in a linear manner. The linear interpolation method is used to calculate the correction factor for each unknown sample.

Delta Calculation: After checking "Mean Correction" or "Interpolation Correction," you can also decide whether to calculate the delta value.



Figure 3.14 External correction strategy

3.3.4 "Export" module

You can select the data to be exported in the "Export" module, which uses the parameters edited in "Formula." Select "Add all" and all of the data ticked in the "Mean" column of "Formula" will be listed. These data can be output to "Workspace (data storage path)." "Export Name" is the final export name, which you can edit again.

Exported Name	Measured Data	Order	
85Rb	85Rb	6	Add
88Sr	88Sr	7	Delete
f-Rb/f-Sr	f-Rb/f-Sr		Add All
87Sr/86Sr	87Sr/86Sr		Delete All
Inner-2SE	Inner-2SE		New Template
Identity	Identity		
K	K		
87Sr/86Sr	Normalized Expr	1	

Figure 3.15 Data export

Iso-Compass provides a formatted report "New Template." Prior to formatting the output, you need to perform the following:

(a) Edit and determine the final "Export Name."

(b) Set the data to be output and the data order in the "Order" column. As indicated in the figure below, the data to be output, including "87Sr/86Sr," "2SE," "84Sr/86Sr," and "Rb/Sr," are numbered in the order of 1, 2, 3...

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)

Exported Name	Measured Data	Order
85Rb	85Rb	6
88Sr	88Sr	7
f-Rb/f-Sr	f-Rb/f-Sr	
87Sr/86Sr	87Sr/86Sr	
Inner-2SE	Inner-2SE	
Identity	Identity	
K	K	
87Sr/86Sr	Normalized Expr	1
2SE	2SE	2
Reference	Reference	
RE(ppm)	RE(ppm)	
84Sr/86Sr	84Sr/86Sr	
Inner-2SE	Inner-2SE	
Identity	Identity	
K	K	
84Sr/86Sr	Normalized Expr	3
2SE	2SE	4
Reference	Reference	
RE(ppm)	RE(ppm)	
83Kr	83Kr	
Rb/Sr	Rb/Sr	5

Figure 3.16 Set data export format and edit data export order

(c) Click "New Template" and set a name for the template.




Figure 3.17 Create new format template

(d) An excel file will be opened automatically. Different sheets have been set for the file, including "Front Cover," "Back Cover," and "Data." You can design the document formats in different sheets according to your needs.

(e) Close the file after saving the Excel file. The template is created successfully and it is saved in the "IsoCompass\template" folder. The data reporting section that follows describes how to use this template.

3.4 Data reduction

After correctly setting the folder path, calibration method template, and sample list of the sample to be processed, click "OK" on the "Create an Analysis Session" page to enter the data reduction page.

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)




Figure 3.18 Data reduction interface

Tool—Option (Fig. 3.19).

The data reduction page includes three major modules:

(a) **Sample List.** The File Name of the Sample List is the name of all “.exp” files in the folder set by the user, and the Sample Name comes from the List file provided by the user. If the name of the “.exp” file does not correspond to the file name in the list, the Sample Name of the “.exp” file will be replaced with the File Name. You can select the file to be processed by ticking the small box in front of the file. Ticking the small box before File Name selects all data.

(b) After selecting a file, the two graphs on the upper right side display the signal distribution of the original sample data and all of the data written during the formula editing. In the original signal map, the two shaded areas are the background signal interval and data signal interval. You can directly select the signal interval manually by dragging the shadow or by entering the required signal interval below. One of the small boxes "All Sync" refers to whether global signal interval selection is performed.

The Formula chart includes a data filtering function. "X" times filtering with "x" sigma means the data culling operation is performed x times. In general, you can tick "All Sync" at the back and fill in "2" times filtering with "2" sigma, which means a 2-times 2x sigma discrete data deletion will be performed on the global sample.

(c) The formula data of all the selected samples are displayed in the lower right part to observe the global data.

You can also enter the Formula module from the data reduction page to make real-time changes to existing calibration methods: Tool—Option.

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)



Figure 3.19 Data reduction interface: Tool toolbar

You can click the "Result Monitor" window in "Analysis Result Watch" in the Tool toolbar to view the correction data in real time. Only data with "Mean," "2SE," and "Count" ticked in the "Formula" module will be displayed in the "Result Monitor" window. The "Result Monitor" window also provides a "Remark" column to record problems with this data during processing.

File Name	Sample	h	Remark	85Rb	88Sr	F-Rb/F-Sr	87Sr/86Sr	2SE	87Sr/86Sr-identity	8i	8i	87Sr/86Sr	8i	84Sr/86Sr	2SE	84Sr/86Sr-identity	B-	B-	84Sr/86Sr	B-	83Kr	Rb/Sr
APR298001	BHVO-2G			0.06486332	2.1440955	0.3231608	0.7033682	0.0001183249	STD			0.703469		-0.000516201	7.195127E-05	STD			0.05657	2.496731E-05	0.03471753	
APR298002	BHVO-2G			0.06504866	2.176949	0.3121229	0.7033635	0.0001424869	STD			0.703469		-0.0003789945	5.837175E-05	STD			0.05657	3.113275E-05	0.03400284	
APR298003	MADA			-2.355541E-06	19.95479		0.7115045	2.452622E-05				0.71187		-0.0007795038	7.984404E-06				0.05657	0.000363019	1.772265E-07	
APR298004	MADA			-2.828848E-05	14.74262		0.7113209	2.533096E-05				0.71187		-0.0011535	7.277403E-06				0.05657	0.000243608	-2.235785E-06	
APR298005	CPX-19			1.986915E-05	0.1508927		0.7007066	0.0009414083						-0.002384551	0.0009283348					5.342895E-06	0.000151381	
APR298006	CPX-20			1.269785E-05	0.126803		0.7029727	0.001244867						-0.003509579	0.00114428					4.336105E-06	0.00011118	
APR298007	CPX-21			4.274479E-05	0.223751		0.6970199	0.0009104258						-0.002402878	0.000558824					7.097415E-05	0.0002132001	
APR298008	CPX-22			1.382162E-06	0.2914331		0.7023777	0.0005360674						-0.002733048	0.000533481					3.515503E-05	4.994646E-06	
APR298009	BHVO-2G			0.08779617	2.995731	0.4416304	0.7032989	0.0001028384	STD			0.703469		-0.0004119181	4.527106E-05	STD			0.05657	7.472655E-05	0.0333556	
APR298010	BHVO-2G			0.08853992	3.075718	0.4305672	0.7032927	9.220439E-05	STD			0.703469		-0.0004194441	4.016672E-05	STD			0.05657	8.218033E-05	0.03295864	
APR298011	BHVO-2G			0.09245032	3.173063	-0.1874822	0.7016424	0.0002118653	STD			0.703469		-0.0003249035	4.502598E-05	STD			0.05657	0.0001483718	0.03333963	
APR298012	BHVO-2G			0.09514986	3.258861	0.5108176	0.7032414	0.0001047585	STD			0.703469		-0.0003176921	4.089286E-05	STD			0.05657	9.693176E-05	0.03300579	
APR298013	BHVO-2G			0.09418997	3.294794	0.3980869	0.7029581	9.47917E-05	STD			0.703469		-0.0003099589	4.262365E-05	STD			0.05657	0.0001117444	0.03236344	
APR298014	BHVO-2G			0.09337777	3.234706	0.4755601	0.7031345	8.778676E-05	STD			0.703469		-0.0003316534	4.891021E-05	STD			0.05657	0.0001021852	0.03304981	
APR298015	BHVO-2G			0.09355755	3.243219	0.5078518	0.7032281	9.070105E-05	STD			0.703469		-0.0004053525	4.354357E-05	STD			0.05657	9.513181E-05	0.03303544	
APR298016	CPX-23			-9.328214E-07	0.3431789		0.7029151	0.0004473072						-0.002363504	0.0003756967					4.0654E-05	-2.812726E-06	
APR298017	CPX-24			-6.436378E-06	0.3737927		0.702795	0.0004544759						-0.002107679	0.0004483183					5.738718E-05	-2.010463E-05	
APR298018	CPX-25			1.104159E-05	0.3739302		0.7027369	0.000444851						-0.001423376	0.0004155483					4.761478E-05	3.262663E-05	
APR298019	CPX-26			2.562385E-05	0.3822149		0.7034827	0.0004287332						-0.001191474	0.0004083392					4.455164E-05	1.700411E-05	
APR298020	CPX-27			3.760124E-05	0.1444246		0.7034919	0.001081239						0.0004596535	0.0008981477					2.408598E-05	0.0002992585	
APR298021	CPX-28			0.0003219305	0.166399		0.7022896	0.001012116						-0.0009495486	0.0007123657					-1.051029E-06	0.00227811	
APR298022	BHVO-2G			0.08459489	3.018714	0.463499	0.7033094	8.86388E-05	STD			0.703469		-0.0004580266	4.349639E-05	STD			0.05657	5.416639E-05	0.03209191	
APR298023	BHVO-2G			0.08154732	3.841909	0.4737147	0.7033133	9.695248E-05	STD			0.703469		-0.000463477	4.906535E-05	STD			0.05657	5.319213E-05	0.03284199	

Figure 3.20 Data display interface after processing

The Export toolbar includes two data export modes: one is the ordinary "Export," which exports all of the data added in the "Export" module, while the other is the templated data export mode, "Formatted Export," Which uses an edited template to output the data. All exported data is saved in "Workspace (data save path)."

3.5 Saving data and exporting reports

All corrected parameters can be saved in "Save Session" in the File toolbar. The saved file can be read directly. However, the saved session can only be observed and cannot be edited. The file is saved in "Workspace (data save path)."

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)




Figure 3.21 Saved session




Figure 3.22 Formatted data export mode

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)




Figure 3.23 Formatted data export cover




Figure 3.24 Formatted data export: data page

4. Laser MC-ICP-MS data reduction

Example: zircon Hf isotope

The zircon Hf isotope is an important research object in micro-area isotope analysis. In this section, the zircon Hf isotope is used as an example to explain to the reader in detail how to use Iso-Compass to carry out micro-area isotope data reduction.

4.1 Data preparation

Prepare a zircon Hf isotope data folder, which contains raw data (".exp" format). All of the isotope signals to be tested are required to appear in the raw file. An edited list of data names in the form of APR27C_LIST.xls is required.




Figure 4.1 Raw data folder

4.2 Data import

Open the Iso-Compass software. Set the working folder path "Workspace," raw data folder path "Folder," file type "File type," and sample name list path "SampleList."

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)




Figure 4.2 Data import

Select the Hf isotope data reduction method "Approach" – "CUG-LA-Hf." After selecting the method, click "Preview/New" to check the method.




Figure 4.3 Checking data reduction

4.3 Data reduction

After confirming that the data reduction method is correct, click "OK" to enter the data reduction page.

On the data reduction page, you need to check the signal interval of each

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)

data point individually to verify the precision of the data and accuracy of the reference substance.




Figure 4.3 Raw data reduction interface

4.4 Data export

After completing the data check, select "Export" – "Formatted Export" and use the established data export template to export the currently processed data in a formatting scheme.




Figure 4.4 Data export

You can then obtain the final formatted data:

Iso-Compass Isotope Data Reduction Software (Ver. 1.0)

中国地质大学 地质过程与矿产资源国家重点实验室 State Key Laboratory of Geological Processes and Mineral Resources		<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>File Name</th> <th>Sample Name</th> <th>Remarks</th> <th>176Hf/177Hf</th> <th>2SE</th> <th>Reference</th> <th>RE(ppm)</th> </tr> </thead> <tbody> <tr><td>APR27C001</td><td>91500</td><td></td><td>0.282295</td><td>0.000030</td><td>0.282302</td><td>-25</td></tr> <tr><td>APR27C002</td><td>91500</td><td></td><td>0.282343</td><td>0.000037</td><td>0.282302</td><td>144</td></tr> <tr><td>APR27C003</td><td>GJ-1</td><td></td><td>0.282041</td><td>0.000024</td><td>0.282013</td><td>101</td></tr> <tr><td>APR27C004</td><td>GJ-1</td><td></td><td>0.282031</td><td>0.000023</td><td>0.282013</td><td>63</td></tr> <tr><td>APR27C005</td><td>TEM</td><td></td><td>0.282684</td><td>0.000024</td><td>0.282686</td><td>-7</td></tr> <tr><td>APR27C006</td><td>TEM</td><td></td><td>0.282683</td><td>0.000027</td><td>0.282686</td><td>-12</td></tr> <tr><td>APR27C007</td><td>S-01</td><td></td><td>0.282478</td><td>0.000036</td><td></td><td></td></tr> <tr><td>APR27C008</td><td>S-02</td><td></td><td>0.282449</td><td>0.000037</td><td></td><td></td></tr> <tr><td>APR27C009</td><td>S-03</td><td></td><td>0.282533</td><td>0.000039</td><td></td><td></td></tr> <tr><td>APR27C010</td><td>S-04</td><td></td><td>0.282454</td><td>0.000045</td><td></td><td></td></tr> <tr><td>APR27C011</td><td>S-05</td><td></td><td>0.282489</td><td>0.000047</td><td></td><td></td></tr> <tr><td>APR27C012</td><td>S-06</td><td></td><td>0.282416</td><td>0.000039</td><td></td><td></td></tr> <tr><td>APR27C013</td><td>S-07</td><td></td><td>0.282493</td><td>0.000047</td><td></td><td></td></tr> <tr><td>APR27C014</td><td>S-08</td><td></td><td>0.282476</td><td>0.000038</td><td></td><td></td></tr> <tr><td>APR27C015</td><td>S-09</td><td></td><td>0.282443</td><td>0.000024</td><td></td><td></td></tr> <tr><td>APR27C016</td><td>S-10</td><td></td><td>0.282448</td><td>0.000042</td><td></td><td></td></tr> <tr><td>APR27C017</td><td>91500</td><td></td><td>0.282286</td><td>0.000024</td><td>0.282302</td><td>-58</td></tr> <tr><td>APR27C018</td><td>91500</td><td></td><td>0.282337</td><td>0.000031</td><td>0.282302</td><td>124</td></tr> <tr><td>APR27C019</td><td>S-11</td><td></td><td>0.282514</td><td>0.000037</td><td></td><td></td></tr> <tr><td>APR27C020</td><td>S-12</td><td></td><td>0.282489</td><td>0.000042</td><td></td><td></td></tr> <tr><td>APR27C021</td><td>S-13</td><td></td><td>0.282513</td><td>0.000037</td><td></td><td></td></tr> <tr><td>APR27C022</td><td>S-14</td><td></td><td>0.282435</td><td>0.000039</td><td></td><td></td></tr> <tr><td>APR27C023</td><td>S-15</td><td></td><td>0.282499</td><td>0.000042</td><td></td><td></td></tr> <tr><td>APR27C024</td><td>S-16</td><td></td><td>0.282466</td><td>0.000036</td><td></td><td></td></tr> <tr><td>APR27C025</td><td>S-17</td><td></td><td>0.282471</td><td>0.000036</td><td></td><td></td></tr> <tr><td>APR27C026</td><td>S-18</td><td></td><td>0.282506</td><td>0.000044</td><td></td><td></td></tr> <tr><td>APR27C027</td><td>S-19</td><td></td><td>0.282486</td><td>0.000034</td><td></td><td></td></tr> <tr><td>APR27C028</td><td>S-20</td><td></td><td>0.282435</td><td>0.000042</td><td></td><td></td></tr> <tr><td>APR27C029</td><td>91500</td><td></td><td>0.282316</td><td>0.000030</td><td>0.282302</td><td>49</td></tr> <tr><td>APR27C030</td><td>91500</td><td></td><td>0.282294</td><td>0.000030</td><td>0.282302</td><td>-29</td></tr> <tr><td>APR27C031</td><td>S-21</td><td></td><td>0.282471</td><td>0.000032</td><td></td><td></td></tr> <tr><td>APR27C032</td><td>S-22</td><td></td><td>0.282483</td><td>0.000035</td><td></td><td></td></tr> <tr><td>APR27C033</td><td>S-23</td><td></td><td>0.282494</td><td>0.000031</td><td></td><td></td></tr> <tr><td>APR27C034</td><td>S-24</td><td></td><td>0.282467</td><td>0.000037</td><td></td><td></td></tr> <tr><td>APR27C035</td><td>S-25</td><td></td><td>0.282470</td><td>0.000037</td><td></td><td></td></tr> <tr><td>APR27C036</td><td>S-26</td><td></td><td>0.282484</td><td>0.000043</td><td></td><td></td></tr> <tr><td>APR27C037</td><td>S-27</td><td></td><td>0.282489</td><td>0.000042</td><td></td><td></td></tr> </tbody> </table>						File Name	Sample Name	Remarks	176Hf/177Hf	2SE	Reference	RE(ppm)	APR27C001	91500		0.282295	0.000030	0.282302	-25	APR27C002	91500		0.282343	0.000037	0.282302	144	APR27C003	GJ-1		0.282041	0.000024	0.282013	101	APR27C004	GJ-1		0.282031	0.000023	0.282013	63	APR27C005	TEM		0.282684	0.000024	0.282686	-7	APR27C006	TEM		0.282683	0.000027	0.282686	-12	APR27C007	S-01		0.282478	0.000036			APR27C008	S-02		0.282449	0.000037			APR27C009	S-03		0.282533	0.000039			APR27C010	S-04		0.282454	0.000045			APR27C011	S-05		0.282489	0.000047			APR27C012	S-06		0.282416	0.000039			APR27C013	S-07		0.282493	0.000047			APR27C014	S-08		0.282476	0.000038			APR27C015	S-09		0.282443	0.000024			APR27C016	S-10		0.282448	0.000042			APR27C017	91500		0.282286	0.000024	0.282302	-58	APR27C018	91500		0.282337	0.000031	0.282302	124	APR27C019	S-11		0.282514	0.000037			APR27C020	S-12		0.282489	0.000042			APR27C021	S-13		0.282513	0.000037			APR27C022	S-14		0.282435	0.000039			APR27C023	S-15		0.282499	0.000042			APR27C024	S-16		0.282466	0.000036			APR27C025	S-17		0.282471	0.000036			APR27C026	S-18		0.282506	0.000044			APR27C027	S-19		0.282486	0.000034			APR27C028	S-20		0.282435	0.000042			APR27C029	91500		0.282316	0.000030	0.282302	49	APR27C030	91500		0.282294	0.000030	0.282302	-29	APR27C031	S-21		0.282471	0.000032			APR27C032	S-22		0.282483	0.000035			APR27C033	S-23		0.282494	0.000031			APR27C034	S-24		0.282467	0.000037			APR27C035	S-25		0.282470	0.000037			APR27C036	S-26		0.282484	0.000043			APR27C037	S-27		0.282489	0.000042		
File Name	Sample Name	Remarks	176Hf/177Hf	2SE	Reference	RE(ppm)																																																																																																																																																																																																																																																																											
APR27C001	91500		0.282295	0.000030	0.282302	-25																																																																																																																																																																																																																																																																											
APR27C002	91500		0.282343	0.000037	0.282302	144																																																																																																																																																																																																																																																																											
APR27C003	GJ-1		0.282041	0.000024	0.282013	101																																																																																																																																																																																																																																																																											
APR27C004	GJ-1		0.282031	0.000023	0.282013	63																																																																																																																																																																																																																																																																											
APR27C005	TEM		0.282684	0.000024	0.282686	-7																																																																																																																																																																																																																																																																											
APR27C006	TEM		0.282683	0.000027	0.282686	-12																																																																																																																																																																																																																																																																											
APR27C007	S-01		0.282478	0.000036																																																																																																																																																																																																																																																																													
APR27C008	S-02		0.282449	0.000037																																																																																																																																																																																																																																																																													
APR27C009	S-03		0.282533	0.000039																																																																																																																																																																																																																																																																													
APR27C010	S-04		0.282454	0.000045																																																																																																																																																																																																																																																																													
APR27C011	S-05		0.282489	0.000047																																																																																																																																																																																																																																																																													
APR27C012	S-06		0.282416	0.000039																																																																																																																																																																																																																																																																													
APR27C013	S-07		0.282493	0.000047																																																																																																																																																																																																																																																																													
APR27C014	S-08		0.282476	0.000038																																																																																																																																																																																																																																																																													
APR27C015	S-09		0.282443	0.000024																																																																																																																																																																																																																																																																													
APR27C016	S-10		0.282448	0.000042																																																																																																																																																																																																																																																																													
APR27C017	91500		0.282286	0.000024	0.282302	-58																																																																																																																																																																																																																																																																											
APR27C018	91500		0.282337	0.000031	0.282302	124																																																																																																																																																																																																																																																																											
APR27C019	S-11		0.282514	0.000037																																																																																																																																																																																																																																																																													
APR27C020	S-12		0.282489	0.000042																																																																																																																																																																																																																																																																													
APR27C021	S-13		0.282513	0.000037																																																																																																																																																																																																																																																																													
APR27C022	S-14		0.282435	0.000039																																																																																																																																																																																																																																																																													
APR27C023	S-15		0.282499	0.000042																																																																																																																																																																																																																																																																													
APR27C024	S-16		0.282466	0.000036																																																																																																																																																																																																																																																																													
APR27C025	S-17		0.282471	0.000036																																																																																																																																																																																																																																																																													
APR27C026	S-18		0.282506	0.000044																																																																																																																																																																																																																																																																													
APR27C027	S-19		0.282486	0.000034																																																																																																																																																																																																																																																																													
APR27C028	S-20		0.282435	0.000042																																																																																																																																																																																																																																																																													
APR27C029	91500		0.282316	0.000030	0.282302	49																																																																																																																																																																																																																																																																											
APR27C030	91500		0.282294	0.000030	0.282302	-29																																																																																																																																																																																																																																																																											
APR27C031	S-21		0.282471	0.000032																																																																																																																																																																																																																																																																													
APR27C032	S-22		0.282483	0.000035																																																																																																																																																																																																																																																																													
APR27C033	S-23		0.282494	0.000031																																																																																																																																																																																																																																																																													
APR27C034	S-24		0.282467	0.000037																																																																																																																																																																																																																																																																													
APR27C035	S-25		0.282470	0.000037																																																																																																																																																																																																																																																																													
APR27C036	S-26		0.282484	0.000043																																																																																																																																																																																																																																																																													
APR27C037	S-27		0.282489	0.000042																																																																																																																																																																																																																																																																													
<h2 style="margin: 0;">检测报告</h2> <p style="margin: 0;">DETECTING AND ANALYZING REPORT</p>																																																																																																																																																																																																																																																																																	
报告编号 SERIES NUMBER 样品名称 SAMPLE NAME 委托单位 ENTRUSTING UNIT 检测类别 DETECTING TYPE 报告制作人 RESPONSIBLE PERSON	MC19_024_ZHANGWEN 磁石 中国地质大学 (武汉) 南区磁石中院位置 张文																																																																																																																																																																																																																																																																																
检测单位 (公章) DETECTING AND ANALYZING UNIT	中国地质大学地质过程与矿产资源国家重点 实验室																																																																																																																																																																																																																																																																																
报告发送日期: REPORTING DATE	2019/8/16																																																																																																																																																																																																																																																																																

Figure 4.5 Formatted data export results