

Attribute	Explanation
VolumeName	The file name.
Manufacturer	Manufacturer of the equipment that produced the Composite Instances.
SeriesDescription	Description of the Series.
ManufacturerModelName	Manufacturer's model name of the equipment that is to be used for beam delivery.
PatientSex	Sex of the named Patient.
PatientAge	Age of the Patient.
ReconstructionDiameter	Diameter, in mm, of the region from within which the data was used in creating the reconstruction of the image. Data may exist outside this region and portions of the patient may exist outside this region.
DistanceSourceToDetector	Distance in mm from source to detector center.
DistanceSourceToPatient	Distance in mm from source to the table, support or bucky side that is closest to the Imaging Subject, as measured along the central ray of the X-Ray beam.
GantryDetectorTilt	Nominal angle of tilt in degrees of the scanning gantry. Not intended for mathematical computations.
TableHeight	The height of the patient table in mm. The range and values of this element are determined by the manufacturer.
RotationDirection	Direction of rotation of the source when relevant, about nearest principal axis of equipment.
ExposureTime	Duration of X-Ray exposure in msec.
XRayTubeCurrent	X-Ray Tube Current in mA.
Exposure	The exposure expressed in mAs, for example calculated from Exposure Time and X-Ray Tube Current.
FilterType	Type of filter(s) inserted into the X-Ray beam (e.g., wedges).
GeneratorPower	Power in kW to the X-Ray generator.
FocalSpots	Nominal focal spot size in mm used to acquire this image.
ConvolutionKernel	A label describing the convolution kernel or algorithm used to reconstruct the data. A single value shall be present.
PatientPosition	Patient position descriptor relative to the equipment.
RevolutionTime	The time in seconds of a complete revolution of the source around the gantry orbit.
SingleCollimationWidth	Adjacent physical detector rows may have been combined to form a single effective acquisition row.
TotalCollimationWidth	This will be equal the number of effective detector rows multiplied by single collimation width.
TableSpeed	The distance in mm that the table moves in one second during the gathering of data that resulted in this image.
TableFeedPerRotation	Motion of the table (in mm) during a complete revolution of the source around the gantry orbit.
SpiralPitchFactor	Ratio of the Table Feed per Rotation to the Total Collimation Width.
DataCollectionCenterPatient	The x, y, and z coordinates (in the Patient-Based Coordinate System) in mm of the center of the region in which data were collected.
ReconstructionTargetCenterPatient	If the reconstructed image is not magnified or panned the value corresponds with the Data Collection Center Attribute.
ExposureModulationType	A multi-valued label describing the type of current modulation used for the purpose of limiting the dose.
CTDIvol	It describes the average dose for this image for the selected CT conditions of operation.
ImagePositionPatient	Image Position specifies the x, y, and z coordinates of the upper left hand corner of the image; it is the center of the first voxel transmitted.
ImageOrientationPatient	Image Orientation specifies the direction cosines of the first row and the first column with respect to the patient.
SliceLocation	Slice Location is defined as the relative position of the image plane expressed in mm. This information is relative to an unspecified implementation specific reference point.
SamplesPerPixel	Samples per Pixel is the number of separate planes in this image.
PhotometricInterpretation	The value of Photometric Interpretation specifies the intended interpretation of the image pixel data.
Rows	Number of rows in the image.
Columns	Number of columns in the image.
XYSpacing	Physical distance between the center of each pixel along the x and y axes, respectively, in the patient.
RescaleIntercept	The value b in relationship between stored values (SV) and the output units.
RescaleSlope	m in the equation specified by Rescale Intercept.
RescaleType	Specifies the output units of Rescale Slope and Rescale Intercept.
NumberofSlices	Number of slices.
ZSpacing	Physical distance between the center of each pixel along the z axis in the patient.
StudyDate	Date the Study started.