

Review

Nanotheranostics

2024; 8(4): 506-520. doi: 10.7150/ntno.96286

Therapeutic and Contrast Agents for Photoacoustic Imaging-Guided Photothermal Therapy: A Narrative Review

Donni Kis Apriyanto^{1,2}, Mitrayana¹, Andreas Setiawan³, Rini Widyaningrum⁴

1. Department of Physics, Faculty of Mathematics and Natural Sciences, Universitas Gadjah Mada, Yogyakarta 55281, Indonesia.

2. Department of Physics, Faculty of Mathematics and Natural Sciences, Universitas Lampung, Bandar Lampung 35141, Indonesia.

3. Department of Physics, Faculty of Science and Mathematics, Universitas Kristen Satya Wacana, Salatiga 50711, Indonesia.

4. Department of Dentomaxillofacial Radiology, Faculty of Dentistry, Universitas Gadjah Mada, Yogyakarta 55281, Indonesia.

🖂 Corresponding author: Rini Widyaningrum, Email: rinihapsara@ugm.ac.id.

© The author(s). This is an open access article distributed under the terms of the Creative Commons Attribution License (https://creativecommons.org/licenses/by/4.0/). See https://ivyspring.com/terms for full terms and conditions.

Received: 2024.03.14; Accepted: 2024.06.29; Published: 2024.08.01

Abstract

Photoacoustic imaging is a hybrid modality that combines high-contrast and spectroscopy-based optical imaging specificity with the high spatial resolution of ultrasonography. This review highlights the development and progress of photoacoustic imaging technology over the past decade. This imaging technology has evolved to be more user-friendly, cost-effective, and portable, demonstrating its potential for diverse clinical applications. A potential clinical application lies in the use of photoacoustic imaging as a guiding tool for photothermal therapy. This review was conducted by initially filtering through three databases, namely, Google Scholar, PubMed, and Scopus, resulting in 460 articles published between 2019 and May 2023. Of these, 54 articles were deemed suitable for review after identification. The selected articles were research papers focusing on the development of therapeutic agents that enhance contrast in photoacoustic imaging. All reviewed articles tested these agents both *in vitro* and *in vivo*. This review focuses on wavelength absorption and radiation sources for photothermal therapy. The developed agents predominantly used NIR-I wavelengths, whereas the NIR-II region has been less explored, indicating significant potential for future research. This review provides comprehensive insights into the advancement of compounds serving as therapeutic agents and contrast agents in photoacoustic imaging-guided photothermal therapy.

Keywords: near-infrared wavelength, radiation, absorption, in vitro, in vivo

Introduction

Photoacoustic imaging is a hybrid modality that combines high-contrast and spectroscopy-based optical imaging specificity with the high spatial resolution of ultrasonography [1]. Several instrument configurations in photoacoustic imaging are employed: optical resolution-photoacoustic microscopy, acoustic resolution-photoacoustic microscopy, photoacoustic tomography/computed tomography, and photoacoustic endoscopy [2]. Each configuration has different resolutions and depths and can be used as required. A review of the development and progress of photoacoustic imaging technology over the past decade revealed that this imaging technology

has evolved to be more user-friendly, cost-effective, and portable, demonstrating its potential for diverse clinical applications [3]. Research involving the application of photoacoustic imaging techniques as monitoring and detection tools in the medical field has been extensive. These include monitoring blood clots during microbubble-assisted sonothrombolysis [4], studying hemodynamic changes caused by blood pressure variations [5], assessing the efficacy of the vascular disrupting agent 5,6-dimethylxanthenone-4-acetic acid through quantitative analysis [6], and studying deep vascularization for regenerative applications [7]. Furthermore, photoacoustic imaging holds potential as a guiding tool for various medical procedures requiring intervention, such as drug delivery, therapy, surgery, and biopsy [8].

applications, In addition to diagnostic photoacoustic research extensively focuses on the development of diagnostic and therapeutic systems (theranostics) and drug monitoring strategies [9]. Photothermal therapy (PTT) is currently under development and uses photoacoustic imaging as a guiding tool. In addition to conventional cancer therapies, PTT represents an alternative and innovative approach to cancer treatment. PTT is a therapeutic modality that uses light-absorbing agents to heat and selectively destroy cancer cells or other abnormal tissues without harming healthy cells [10]. The hybridization of photoacoustic imaging and PTT demonstrates clear advantages in monitoring and treating tumors because of its high efficiency, spatial specificity, fewer side effects, and ease of operation [11].

Photoacoustic imaging has been extensively investigated as a guide for PTT. These studies have focused on various compounds that assist in disease therapy and agents to enhance the contrast of photoacoustic imaging results. These agents assist in the more precise visualization of targets such as tumors or other diseases, especially at the cellular Using photoacoustic imaging, medical level. professionals can obtain valuable information about tumors' location, size, and characteristics. In addition to their use in diagnosis, these agents possess therapeutic properties. They can generate heat when exposed to light radiation, which can be utilized to destroy or damage cancer cells in photothermal therapy. This approach allows for more targeted and precise treatment, minimizing damage to surrounding healthy tissue. Compounds used as agents in PTT are formulated using various compound mixtures; as a result, they possess different maximum absorption wavelengths. Therefore, this review summarizes recent advancements in agents for photoacoustic imaging as guides for PTT based on their components, absorption wavelength ranges, and in vitro and in vivo PTT applications.

Furthermore, this review can serve as a reference for the development of agents that enhance image contrast in photoacoustic imaging and possess therapeutic properties for treating diseases such as cancer. Thus, these potential agents can be highly efficient and effective diagnostic and therapeutic tools. This review also provides conclusions and a brief discussion on the challenges and perspectives and the potential future directions of PTT agents that simultaneously enhance photoacoustic imaging contrast.

Methods

А comprehensive literature search was conducted to filter articles related to photoacoustic imaging-guided PTT. Articles were obtained from the Google Scholar, PubMed, and Scopus databases. The literature search used the keywords "photoacoustic imaging" and "photothermal therapy." The Boolean operator "AND" and phrase searching (" ") were explicitly used to facilitate the search. Clinical trials, case series, case reports, and other relevant articles were considered for review. The selected articles were those published between early 2019 and May 2023. The identified review articles should have been included in the initial evaluation. The initial assessment involved identifying titles through literature searches. Articles not written in English or inaccessible in full text were excluded from the initial evaluation. A flowchart of the search and literature selection process is shown in Figure 1. The included studies were evaluated and then organized in a table based on the article title, agent and compound type, cells used in vitro or in vivo, peak absorption wavelength, photoacoustic wavelength, PTT laser wavelength, and compound temperature change when exposed to laser irradiation.

The exclusion criteria were applied after creating the table to eliminate articles. Further selection of articles was based on titles containing "photoacoustic imaging," "guided," and "photothermal therapy." In addition, articles featuring compounds used as therapeutic and contrast agents were considered. Articles combining *in vitro* and *in vivo* methods were selected for inclusion in this literature review.

Results and Discussion

Initial screening across the three databases vielded 460 articles. After removing duplicates, 345 articles were included in the initial evaluation. Based on further selection using the following inclusion criteria, 95 articles were selected for further review. During the subsequent review process, 54 articles were deemed suitable for review (Figure 1). The summary of the article review is organized into Table 1, Table 2, Table 3, and Table 4. Table 1 presents the optical characteristics of the agent compounds, including the agent name, type, absorption wavelength (λ_{abs}), wavelength used for generating photoacoustic images (λ_{PA}), imaging modality used, and tool used for image generation. Table 2 summarizes the distribution of agents based on their absorption wavelengths. Table 3 summarizes the thermal characteristics of the agent compounds in the in vitro methods. This includes the agent's name, cells used in vitro, therapeutic technique employed,

radiation wavelength delivered (λ_{Rad}), radiation source power, duration of agent radiation, temperature increase, and temperature achieved during the provided radiation time, and photothermal conversion efficiency of the agent. **Table 4** presents data related to the thermal characteristics of the agent compounds in the *in vivo* methods. It contains the agent's name, cells implanted in experimental animals, therapeutic technique used, tool for visualizing the agent's temperature in experimental animals, radiation wavelength delivered (λ_{Rad}), radiation source power, duration of agent radiation, temperature increase, and temperature achieved during the provided radiation time.

Optical characteristics of the agents

This review categorizes agents into organic, inorganic, hybrid organic-inorganic, organic dye, organic-inorganic dye, and dye categories. Organic or inorganic refers to compounds formed on basis of the organic or inorganic components, whereas hybrid organic-inorganic refers to mixed compounds formed from organic and inorganic components. Organic and organic-inorganic dyes are colorants derived from organic components or a mixture of organic-inorganic components. This classification is based on the information in the articles or material component of the agents used in the studies reviewed in this work.

Figure 2 indicates that the compounds used in photoacoustic imaging-guided PTT are mainly

derived from a combination of organic-inorganic components, with compounds originating from organic and inorganic elements following closely. The prevalent use of organic materials stems from their potential for biocompatibility and biodegradability within biological systems compared with inorganic materials. Organic materials also exhibit high absorption in the far-red or near-infrared (NIR) region, presenting significant potential for PTT [12]. As shown in **Table 1** and **Figure 2**, studies employing dyes as contrast and therapeutic agents are fewer than those employing other agent types.

One of the main challenges in photoacoustic imaging is the limited depth of light penetration in biological tissues, especially in the visible and NIR-I Additionally, regions. some non-contact photoacoustic imaging methods may encounter constraints related to signal travel distance, especially when imaging deep tissues. These constraints can lead to lower resolution in deep tissues than direct contact imaging methods (Ying et al., 2024). These limitations reduce the effectiveness of imaging deep tissues and organs, limiting their clinical applications. Therefore, to enhance treatment effectiveness, developing techniques capable of overcoming these barriers is necessary, such as using radiation sources with longer wavelengths or developing contrast agents with optimized optical properties for better tissue penetration.



Figure 1. Selection diagram and literature search.

Agent	Туре	λ_{abs} (nm)	λ_{PA} (nm)	Imaging modality	PAI tools	Ref
Polydopamine-coated mesoporous silica-gold nanorods (PDA-AuNRs@MSN)	Organic-inorganic	450	875	PAI	Vevo LAZR	[15]
Diketopyrrolopyrrole benzothiadiazole nanoparticles (DPP-BT NP)	Organic Semiconducting	686	808; 730	Fluorescent; PA	Nexus 128	[16]
Polyethylene glycol- and polyethyleneimine-modified gold nanorods (mPEG-PEI-AuNRs); mPEG-PEI and calcium carbonate nanonarticles (mPEG-PEI/CaNPs)	Organic-inorganic	n/a	n/a	Fluorescent; PAI	MSOT InVision 128	[17]
Gold, 11-mercaptoundecanoic acid, and 10-mercaptodecyl) trimethylammonium bromide	Organic-inorganic	524	n/a	PAI	MSOT InVision 128	[18]
Benzobisthiadiazole (BBTD-1302)	Organic semiconducting polymer	942	1302; 942	Fluorescent; PAI	MSOT imaging system	[19]
Emeraldine salt of polyaniline in sodium bis(2-ethylhexyl) sulfosuccinate (PANI-ES@AOT)	Organic-inorganic	n/a	808; 970	PAI	Vero LAZR-X mX550D	[20]
1T Molybdenum disulfide (1T-MoS ₂); 2H Molybdenum disulfide (2H-MoS ₂)	Inorganic	1300; 427, 606, 653	1280	PAI	n/a	[21]
Benzodithiophene-Thieno[3,4-b]thiophene (BDT-TH); Benzodithiophene-Naphthalenediimide (BDT-NDI)	Organic Semiconducting polymer	516; 700	532; 700	PAT	Homemade PA system	[22]
Boron-dipyrromethene-Fe nanoparticles (BDP-Fe NPs)	Organic-inorganic	1300	1064	PAI	n/a	[23]
3-azidopropylamine (-N ₃ ⁻) with Bi ₂ S ₃ -Ag ₂ S-DATS@BSA nanoparticles (AB-DS@BSA-N ₃ NYs)	Organic-inorganic	900	808	Fluorescent; PAI	MSOT; II900/1700 in vivo imaging System (808 nm, 150 ms)	[24]
Hollow mesoporous carbon spheres-polyethylene glycol-gambogic acid (HMCS-PEG-GA)	Organic-inorganic	808	715	PAI	High-resolution preclinical photoacoustic system	[25]
Melanin nanoparticles H ₂ (MNPH ₂)	Organic	457	808	PAI	The MSOT scanner consists of a range of 680– 980 nm	[26]
Folic acid-conjugated chitosan-functionalized graphene oxide (FA-CS-GO)	Organic-inorganic	578, 720	532	PAI	PAM system (using a Q-switched diode-pumped solid-state laser (SPOT-10-100-532) and a commercially focused transducer (V324-SM); a reflection-mode fluorescent imaging system	[27]
Diketopyrrolopyrrole (DPP) was copolymerized with strong electron-donating substitutes (fluorene and thiophene) (PFTDPP); PFTDPP combined with S-nitroso-N-acetylpenicillamine (PFTDPP-SNAP NPs)	Organic	350–500, 550–900	808	Fluorescent; PAI	PA imaging	[28]
Porphyrin-based covalent organic framework	Organic	590, 297	n/a	PAI	MSOT InVision 128	[29]
CoCuFe-selenide (CCFS)-polyvinyl pyrrolidone (PVP)-L-arginine (L-Arg) (CPA nanocomposite)	Organic	808	n/a	PAI	MSOT InVision 128	[30]
Squaraine dye nanoparticles (SQ-NPs)	Dye	600-800	760	PAI	n/a	[31]
Molybdenum(V)-based polyoxometalate (POMs)	Inorganic	1065	1064	PAI	MSOT InVision 128	[32]
Cyclo(Arg-Gly-Asp-Dphe-Lys) @ polymer dithienopyrrole -thiadiazolobenzotriazole (cRGD@PT NPs)	Organic	1000	1064	PAI	LOS-3D Photoacoustic imaging	[33]
ISSzyme	Organic-inorganic	700	532	PAI	Homemade PACT system (Q-switched Nd: YAG laser (532 nm))	[34]
SY1080 NPs	Organic dye	820	808	Fluorescent; PAI	MSOT imaging system	[35]
Bodipy-1O2 methoxyl (B-OMe-NPs)	Organic	735	809	Fluorescent; PAI	n/a	[36]
PEG-modified boron quantum dots (BQDs-PEG)	Organic-inorganic	n/a	808	Fluorescent; PAI	n/a	[37]
Y16-Pr-PEG NPs	Organic	800	808	Fluorescent; PAI	Confocal laser scanning microscopy (Leica SP8); Near-infrared InGaAs array imaging detector (EX, 808 nm; 2 W/cm2)	[38]
Gold nanostar polyaniline hyaluronic acid (AuNSPHs) pyrazino[2,3g]quinoxaline (PQs)-TPA NPs	Organic-inorganic Organic small	n/a 675	850 760	PAI Fluorescent; PAI	VisualSonics Vevo LAZR MSOT InVision 128	[39] [40]
Melanin-nanoparticular polydopamine-indocyanine	molecules Organic dye	808	808	Fluorescent; PAI	VEVO LAZR-X (at a wavelength of 680-970 nm	[41]
Pt/Te nanoheterostructures (PT); mPT	Inorganic	n/a	744	PAI	744 nm, 120 nJ laser by the high-resolution PAI	[42]
Hemicyanine dye nanoparticles (M1 NPs)	Organic dye	734,	734;	Fluorescent; PAI	Perkin Elmer IVIS Lumina III; PACT scanner	[43]
DEC vioted Ch 925 NDHa	Inorgania	1040 925	1040 925	Eluorosconti DAI	(IIIVISIOII256-1F)	[44]
Au-MBA denotes 4-mercaptobenzoic acid-Cy7	Organic-inorganic	855 497, 590,	835 808	Fluorescent; FAI	n/a n/a	[44] [45]
(Au ₄₄ MBA ₂₆ -Cy7) NCs	dye	794	(00	PAI		
Janus chitosan/ gold quantum dots (J-Au-CS-QD) Gold nanoparticle-coated reduced graphene oxide (rGO-AuNP)	Organic–inorganic Organic–inorganic	808 1000	690 1250	Fluorescent; PAI US + PAI	n/a VisualSonic Vevo 2100 LAZR (40 MHz)	[46] [47]
MoSe _{2-x} nanoflowers (MNFs)	Inorganic	n/a	800	PAI	Endra Nexus 128	[48]
W ₁₈ O ₄₉ nanorods	Inorganic	n/a	n/a	PAI	MSOT imaging system	[49]
Polymer PIIGDTS nanoparticles with folic acid modification (PD-FA)	Organic	450, 845	845	PAI	MSPAM system (wavelength tunable nanosecond pulsed OPO laser (Continuum, USA), ultrasound transducer with a central frequency of 25 MHz (Olympus, USA))	[50]

Agent	Туре	λ_{abs} (nm)	λ_{PA} (nm)	Imaging modality	PAI tools	Ref
Diketopyrrolopyrrole (DPP)-SO; DPP-SS; DPP-SSe NPs	Organic	900; 912; 960		PAI	MSOT imaging system	[51]
${\rm Fe^{3+}}$ complex Zwitterionic dopamine sulfonate (Fe-ZDS)	Organic-inorganic	450-900	680, 808	MRI; PAI	PA imager under the excitation of 680 or 808-nm pulse lasers	ı [52]
Mo ₂ C-derived polyoxometalate (POM)	Inorganic	1060	n/a	PAI	NIR-II PAI images	[53]
DPP-SS; DPP-OF; DPP-SF; DPP-SeF NP	Organic	726; 735; 816; 820	n/a	PAI	MSOT InVision 128	[54]
PEG-modified polydopamine (MPDA-PEG)	Organic	n/a	n/a	PAI	High-resolution preclinical PAI system	[55]
TMB-F4TCNQ; TMB-TCNQ	Organic dye	400; 392	1300	PAI	PA images under 1300-nm laser excitation	[56]
IDCIC NPs	Organic	682	780; – 1000	Fluorescent; PAI	FLI (FC-W-980-30W); PACT Scanner (Invision256-TF)	[57]
Glypican-3-binding peptide-Fe ₃ O ₄ Core/Au shell nanocomplex (GBP-FANB)	Inorganic	520-680	640	Fluorescent; PAI	FLI (IVIS Lumina II (Excitation Filter: 640 nm, Emission Filter: 700 nm); PA (Endra Nexus128)	[58]
Indacenodithiophene (IDT)-BTzTD Pdots	Organic	800	n/a	CT, fluorescent; PAI	SuperArgus (with the following parameters: 50 kV, 300 µA, 33-µm resolution, and 300 ms of exposure time) PAI: homemade PAI system	[59]
Au@Cu ₂ -xSe NPs	Inorganic	NIR	808	CT; PAI	VEVO LAZR-X	[60]
FA-doxorubicin (DOX)-indocyanine green (ICG)- perfluoropentane (PFP)@Lip	Organic	450	n/a	US; PAI	n/a	[61]
Manganese pentacarbonyl bromide-loaded g-carbon nitride/polypyrrole (MnCO@CNPpy)	Organic-inorganic	410-430	n/a	PAI	Vevo LAZR-X (40 MHz)	[62]
Gold nanostars (GNS); Silica@GNS	Organic-inorganic	850	850	PAI	PA system (OPO laser 850 nm; hydrophone (HNC-1500, ONDA, Sunnyvale, CA); amplifier (5072PR, OLYMPUS, Waltham, MA); transducer (Immersion Transducers, Olympus Corp., Tokyo, Japan) with a center frequency of 5 MHz; TWPA-I	[63]
ZrO _{2-x} @PEG/Ce6 (ZPC NPs)	Organic-inorganic	~405; ~660	n/a	Fluorescent; PAI	Fusion FX Spectra; Vevo LAZR-X multimodal imaging system	[64]
Gold nanorod melanin nanoparticles (GNR@MNPs)-Gd-64 Cu	Organic-inorganic	700-900	680-980	PET; MRI; PAI	Siemens Inveon micro PET; 3.0-T MR scanner; MSOT inSight 128 (under a 680–980 nm pulse laser);	[65]
PFTQ-PEG-Gd NPs	Organic-inorganic	760	760	FLI; MRI; PAI	CDD: NIRvana TE 640; micro-MRI scanner (0.5 T, NIUMAG, NMI20-015 V-I, TE/TR = 18/120 ms); Vevo LAZR PAI System	[66]
Ag ₂ S@M-P-RGD	Organic-inorganic	275	744	FLI; PAI	n/a	[67]
Au@PB@Cu2O@BPQDs/PAH NCs	Organic-inorganic	650-950	514; 650	MRI; PAI	3.0 T clinical MRI system; PAI n/a	[68]

PS: n/s signifies that the information provided is not written or not detailed.

Table 2. Summary of agent distribution based on their absorbance wavelength regions

λ _{abs} (nm)	Agents
Visible region (400-700 nm)	PDA-AuNRs@MSN [15]; DPP-BT NP [16]; Au-MUA-TMA [18]; 2H-MoS ₂ [21]; BDT-TH; BDT-NDI [22]; MNPH ₂ [26]; FA-CS-GO [27];
	PFTDPP [28]; COF-366 NPs [29]; ISSzyme [34]; PQs-TPA NPs [40]; Au44MBA26-Cy7 NCs [45]; PD-FA [50]; Fe-ZDS [52]; TMB-F4TCNQ;
	TMB-TCNQ [56]; IDCIC NPs [57]; GBP-FANB [58]; DOX- ICG-PFP@Lip [61]; MnCO@CNPpy [62]; ZPC NPs [64]; Ag:S@M-P-RGD [67];
NIR-I (700-1000 nm)	BBTD-1302 [19]; AB-DS@BSA-N3 NYs [24]; HMCS-PEG-GA [25]; FA-CS-GO [27]; PFTDPP-SNAP NPs [28]; CPA nanocomposite [30];
	SQ-NPs [31]; SY1080 NPs [35]; B-OMe-NPs [36]; Y16-Pr-PEG NPs [38]; Mel-pDA-ICG [41]; M1 NPs [43]; PEGylated Sb-835 NPHs [44];
	Au44MBA26-Cy7 NCs [45]; J-Au-CS-QD [46]; PD-FA [50]; DPP-SO; DPP-SS; DPP-SSe NPs [51]; Fe-ZDS [52]; DPP-SS; DPP-OF; DPP-SF;
	DPP-SeF NP [54]; IDT-BTzTD Pdots [59]; Au@Cu _{2-x} Se NPs [60]; GNS; Silica@GNS [63]; GNR@MNPs-Gd-64 Cu [65]; PFTQ-PEG-Gd NPs
	[66]; Au@PB@Cu ₂ O@BPQDs/PAH NCs [68]
NIR-II (1000-1700 nm)	1T-MoS₂ [21]; BDP-Fe NPs [23]; POMs [32]; cRGD@PT NPs [33]; M1 NPs [43]; rGO-AuNP [47]; POM [53];



Figure 2. Distribution of contrast agent types, including organic, inorganic, hybrid organic-inorganic, organic dye, organic-inorganic dye, and dye.

The wavelength spectrum regions of NIR can be divided into three channels: NIR-I (700–900 nm), NIR-IIa (1300–1400 nm), and NIR-IIb (1500–1700 nm) [13]. In this review, they were categorized into the visible region (400–700 nm), NIR-I (700–1000 nm), and NIR-II (1000–1700 nm) [14] because of the presence of agents with absorption wavelengths in the visible light region.

Table 2 summarizes the distribution of agents based on their absorption wavelengths, which is illustrated in Figure 3. In Figure 3 presents the development of agents that absorb at specific wavelengths. More agents are designed to be absorbed in the NIR-I wavelength region. However, studies using NIR-II are still fewer than those utilizing other wavelength ranges. The use of the NIR-II window (1000-1700 nm) was limited because of the absence of contrast agents with absorption capabilities in the NIR-II region [69]. They reported synthesizing photoacoustic contrast agents that could be absorbed in both the NIR-I (750 nm) and NIR-II (1064 nm) comparing regions, directly the resulting photoacoustic images. The findings indicated that the signal-to-noise ratio of the photoacoustic images generated at a wavelength of 1064 nm was higher than that at 750 nm.



Figure 3. Distribution of the wavelength absorbance of the agents.

Figure 4 displays the distribution of wavelength usage as a radiation source for generating photoacoustic images. In this figure, the NIR-I region is more extensively used as a radiation source. The number of agents in each absorption region in Figure 3 is mostly the same as that in each wavelength region used to generate photoacoustic images in Figure 4. This indicates that the selection of wavelengths as a source for photoacoustic imaging is adjusted according to the agents' absorption wavelengths to maximize the obtained images. Moreover, the choice of photoacoustic imaging wavelengths is also determined by the depth of the object being imaged.

PTT using in vitro methods

Table 3 summarizes the thermal characteristics of the agent compounds in the *in vitro* methods, including the cells used *in vitro*, laser irradiation, agent temperature increase, and agent photothermal conversion efficiency.





Figure 5 shows that 4T1 and HeLa cells are more frequently used than other cell types. The 4T1 cell line is a highly tumorigenic and invasive transplantable tumor cell line [70]. Unlike many tumor models, 4T1 can spontaneously metastasize from the primary tumor in the mammary gland to several distant sites, including the lymph nodes, blood, liver, lungs, brain, and bones.

In addition to serving as contrast agents, these agents are used as cancer therapies through PTT. Therefore, these agents should reach temperatures sufficient to eliminate cancerous cells. Thermal cancer predominantly identified therapy is with temperatures ranging from 40 to 48°C sustained at the treated site for an hour or longer [71]. Moreover, temperatures above 42°C lead to cell death in various tissues [72]. Agents undergo laser irradiation for a specified duration to attain temperatures that can eradicate cancer cells. We cannot directly compare the effectiveness of the agents in temperature elevation because of varying radiation treatments (wavelengths, power, and duration). Figure 6 illustrates the rate of temperature increase calculated from the temperature change data and radiation duration applied to the agents. The highest rate of temperature increase was achieved by rGO-AuNP at 12.2 °C/min using a radiation wavelength of 1061 nm (0.25 W/cm², 300 s) [47]. The exact reason why this agent exhibits a higher rate of temperature increase than other agents is not explicitly stated. Comparatively, for agents irradiated

at the same radiation wavelength (NIR-II), Mo_2C -derived polyoxometalate (POM) irradiated at 1060 nm (1 W/cm², 300 s) achieves a rate of 7.4°C/min [53], and 2H-MoS₂ irradiated at 1064 nm (1 W/cm², 480 s) achieves a rate of 1.55 °C/min [21]. Among agents irradiated in the NIR-I or visible light region, the highest rate of increase is observed in the PQs-TPA NPs agent irradiated at a wavelength of 808

nm (0.8 W/cm², 300 s) at 9.2 °C/min [40]. In the comparison of agent composition, rGO-AuNP containing Au, AuNSPHs, Au₄₄MBA₂₆-Cy₇ NCs, Au@Cu₂-xSe NPs, J-Au-CS-QD, and Au-MUA-TMA exhibited significantly lower rates at 4°C/min, 6.7°C/min, 4°C/min, 2.83°C/min, and 2.8°C/min [18,39,45,46,60].

Table 3. Summary of the effects of photothermal agents in vitro

Agents	Cells	Therapy	λ_{Rad} (nm)	Power Rad (W/cm ²)	Time Rad (s)	ΔT (°C)	T _{max} (°C)	Efficiency (%)	Ref
PDA-AuNRs@MSN	4T1	PTT	808	1	600	31.4	n/a	34.71	[15]
DPP-BT NP	HeLa	PDT; PTT	730	1	600	21.8	43.6	50	[16]
mPEG-PEI-AuNRs; mPEG-PEI/CaNPs	MCF-7, HeLa, 293T	UST; PTT	808	1.5; 0.8	600; 360	n/a	60	27.7	[17]
Au-MUA-TMA	U87MG	PTT	808	1	600	28	65	n/a	[18]
BBTD-1302	4T1	PTT	980	1	n/a	23	n/a	28.6	[19]
PANI-ES@AOT	HeLa, MCF-7	PTT	1064	1	300	31.1	n/a	43.9	[20]
1T-MoS ₂ nanodots: 2H-MoS ₂	A549: HeLa	PTT	808:1064	1	480	43:12.4	n/a	43.3: 21.3	[21]
BDT-TH: BDT-NDI	MCF-7	(LED) PTT	660	0.2	600	25-30	n/a	40	[22]
BDP-Fe NPs	HeI a	PTT	730	n/a	600	23	n/a	49	[23]
AB-DS@BSA-Na NVs	Hen-2	PTT	808	1	300	n/a	60.9	31.6	[24]
HMCS PEC CA	HopC2	PTT	808	1	300	323	573	47.3	[25]
MNPH2	Hop 2	PTT	808	0.75	300	52.5 n/a	63.8	25.7	[26]
	MDA MR 221	DTT	000	0.75	200	n/a	51.4	23.7	[20]
DETODD, DETODD CNAD ND-	MCE 7	DTT	000	1	200	11/a	51.4	11/ a	[20]
COE 2((NP-	MCF-7	I'II DDT. DTT	000	1	200	n/a	57	40	[20]
CDA		PDI; PII	000	0.2; 1.5	500	23 1 0 0	60 50.1	15.07	[29]
CPA nanocomposite	U87mg, HeLa, HepG2, COS7	PII	808	1	600	12.9- 45.9	50.1	72 (pH 7.4) 81 (pH 5.4)	[30]
SQ-NPs	MCF-7; A549	PIT	760	0.5	600	25.1	n/a	52.8	[31]
Molybdenum(V)-based POMs	411	PIT	1064	0.8	300	n/a	64.6	n/a	[32]
cRGD@PT NPs	U87MG, C6, GL261, HU-VECs, and HeLa	PIT	1064	1	300	20	n/a	43.15	[33]
ISSzyme	4T1; HU-VECs	PTT	808	1	300	32.2	54.17	70.9	[34]
SY1080 NPs	4T1	PTT	808	1	300	n/a	52.7	22.3	[35]
B-OMe-NPs	4T1	PTT	735	1	600	n/a	68	66.5	[36]
BQDs-PEG	4T1	PTT	808	1	600	n/a	52	23.7	[37]
Y16-Pr-PEG NPs	4T1	PTT	808	1	300	25	50	82.4	[38]
AuNSPHs	4T1; bEnd.3	PTT	809	1	300	20	n/a	78.6	[39]
PQs-TPA NPs	L929, HeLa, 4T1	PTT	808	0.8	300	46	73.8	59.5	[40]
Mel-pDA-ICG	4T1, CT26	PTT	808	1.2	180	28	n/a	n/a	[41]
PT; mPT	HU-VECs, 4T1	PTT	808	1.5	600	45	n/a	51.84	[42]
M1 NPs	4T1	PTT	808	1	600	n/a	64	77.5	[43]
PEGylated Sb-835 NPHs	4T1,	PTT	808	1	540	60	n/a	62.1	[44]
Au ₄₄ MBA ₂₆ -Cy7 NCs	HeLa,	PTT	808	1.5	360	40	70	65.12	[45]
J-Au-CS-QD	4T1, HEK293	PTT	808	0.5	600	28.3	n/a	43	[46]
rGO-AuNP	SKOV-3	PTT	1061	0.25	300	61	n/a	n/a	[47]
MNFs	L929; 4T1	PTT	808	1	600	50.9	n/a	41.7	[48]
W ₁₈ O ₄₉ nanorods	4T1; HU-VECs	CDT; TT	n/a	n/a	n/a	16.6	n/a	n/a	[49]
PD-FA nanoparticles	HeLa: cos7	PTT	, 845	1	600	n/a	80	62.6	[50]
DPP-SO: DPP-SS: DPP-SSe NPs	A549; H446	PTT	808	0.8	300	35	n/a	79.3	[51]
Fe-ZDS	EMT-6; MCF-7	PTT	808	1.65	300	40	n/a	30.99	[52]
POM	HUVEC, HeLa	CDT: PTT	1060	1	300	37	n/a	47.8	[53]
DPP-SS; DPP-OF; DPP-SF; DPP-SeF NP	A549; H446	PTT	808	0.75	300	-; -; -; 37	50, 52, 56, 62	62	[54]
MPDA-PEG	HEK-293T: A549	PTT	808	2	300	41.5	66.5	41	[55]
TMB-F4TCNO: TMB-TCNO	4T1, HeLa	PTT	1060	1	300	n/a	n/a	`: 4 8	[56]
IDCIC NPs	4T1	PTT	808	1	300	32.6	59.9	78.9	[57]
GBP-FANB	HepG2: PC3	PTT	630	1	600	n/a	67	n/a	[58]
IDT-BTzTD Pdots	HeLa	PTT	808	0.8	600	32.1	574	53.9	[59]
Au@Cua-xSe NPs	4T1	PTT	808	1	300	20	45	56.6	[60]
FA-DOX-ICG-PEP@I in	Y79. Hu-VECs	PTT	808	1	300	>40	>60	n/a	[61]
MnCO@CNPpy	Hela MCE-7 L02	PTT	808	1	300	40	n/a	34.6	[62]
GNS: Silica@GNS	MSC	PTT	808	- 1	900	55	80	n/a	[63]
ZPC NPc	HeI a: 4T1	 PDT- PTT	660:808	0.76:2	300	n/a	69.1	30.99	[64]
CNR@MNPs-Cd-64 Cu	Hon-2	PTT	808	1	300	n/a	60	45.3	[0±] [65]
PETO-PEC-Cd NPe	/T1	PTT	808	1	600-300	n/a	>60	26	[66]
Ag-S@M-P-RCD	HeI a	PTT	808	2	600	n/a	65	28 35	[67]
Au@PB@CuoO@BPODs/PAH NCs	Hela	PDT· PTT	650	- 15	300	n/a	52 5	25.73	[68]
Inter Decuzoebi QD6/171111005	******			1.0	500		54.0		[00]

PS: n/a signifies that the information provided is not written or not detailed.

Figure 7 illustrates the distribution of the maximum temperature increases for the agents. Referring to the 40°C–48°C temperature range for cancer therapy [71], nearly all these agents can be applied for cancer treatment. Silica@GNS achieved the highest temperature at 80°C when irradiated with a wavelength of 808 nm (1 W/cm^2 , 900 s) [63]. Once again, we could not find a precise reason for this agent reaching a higher temperature than others. This necessitates a review of numerous factors such as pH, density, size, and even the composition of the agents.



Figure 8 shows the thermal conversion efficiency distribution based on the agent type. Photothermal conversion efficiency refers to the ability of an agent to convert light energy into heat. This graph compares conversion efficiency among various materials such as organic, inorganic, hybrid, or dyes. This data is essential for understanding how each component enhances the effectiveness of photothermal applications in cancer therapy or local heating. Figure 9 shows the distribution of thermal conversion

efficiency across different absorbance wavelength regions. This graph is used to determine the optimal efficiency at specific absorbance wavelength regions for medical applications. From these figures, we cannot conclude that the agent's components and the length of the radiation laser source directly affect the increase in photothermal conversion efficiency. Improving photothermal conversion facilitates enhanced efficiency in PTT and contributes to an increase in photoacoustic signals [18]. Figure 10 illustrates the relationship between photothermal conversion efficiency and the temperature change produced by the agent. The graph does not show a linear relationship between higher conversion efficiency and temperature increase. Figure 11 depicts the relationship between photothermal conversion efficiency and the temperature the agent achieves. This graph provides insights into the potential of each agent to reach the high temperatures necessary for



Figure 6. Distribution of temperature increase in vitro for the agents.

therapeutic applications.



Figure 7. Distribution of the achieved temperature increase by the agents.



This information helps evaluate and compare agents' performance under intensive use conditions. Various parameters, such as the laser intensity applied, can influence the efficiency and temperature increase. The intensity of the laser light affects the amount of energy absorbed by the photosensitizer or nanoparticles, which in turn influences the rise in temperature. Therefore, essential factors affecting photothermal performance and the methods and mechanisms of laser-induced photothermal energy must be studied [73].



Additionally, the volume fraction, shape, size, and material type of the agent compounds influence the efficiency and temperature rise achievable by the agent [74]. The highest photothermal conversion efficiency is observed in Y16–Pr–PEG NPs, reaching 82.4%, irradiated at a wavelength of 808 nm (1 W/cm2, 300 s), resulting in a temperature increase of 25°C and achieving a temperature of 50°C [38]. In comparison, MPDA-PEG, irradiated at a wavelength of 808 nm (2 W/cm2, 300 s) [55], exhibits a lower photothermal conversion efficiency of 41%, a temperature increase of 41.5°C, and a temperature of 66.5°C.



Although photothermal therapy is promising as an exact treatment method, challenges remain in improving photothermal conversion efficiency and achieving uniform heat distribution within the tumor. Uneven heat distribution can result in incomplete tumor ablation and potential damage to surrounding healthy tissues. Further research is needed to develop new materials with optimized optical and thermal characteristics to enhance the efficiency of photothermal therapy. These parameters include heating conditions (laser intensity, wavelength, radius, and heating time), nanoparticle parameters (volume fraction, shape, size, and material), and boundary conditions (cooling methods and cooling time), all of which are applied to control therapeutic outcomes (Ren et al., 2022).





Figure 12. Distribution of cells used in vivo.

PTT using in vivo methods

Table 4 summarizes the thermal characteristics of the agent compounds in the *in vivo* methods, including the implanted cells in experimental animals, thermal camera tools, laser irradiation, and temperature increase of the agents.

Figure 12 shows that the 4T1 cell line is more frequently used in experimental animals than other cell lines. The 4T1 tumor exhibits several characteristics that make it a suitable animal model for human breast cancer [70]. The progressive metastatic spread of 4T1 to lymph nodes and other organs is highly reminiscent of human breast cancer. The 4T1 tumor is also resistant to 6-thioguanine. This characteristic enables accurate counting of metastatic cells, even when spread at submicroscopic levels in distant organs.

Similarly, as in *in vitro* studies, agents are used as a therapeutic agent for cancer implanted in experimental animals. These agents are injected into the experimental animals, accumulate in the cancer cells, and are then exposed to radiation, causing a sufficient increase in temperature to eliminate cancer cells. Figure 13 illustrates the distribution of the rate of temperature increase of the agents in vivo, whereas Figure 14 displays the distribution of the achieved temperature of the agents in vivo. Fe-ZDS at 7.5°C/min recorded the highest rate, with radiation at a wavelength of 808 nm (0.82 W/cm², 360 s) [52]. Meanwhile, $1T-MoS_2$ nanodots achieved the maximum temperature increase at 65.7°C using radiation at a wavelength of 808 nm $(1 \text{ W/cm}^2, 480 \text{ s})$. Figure 14 and Table 4 also indicate agents that can achieve a temperature range of 40°C-48°C, which is helpful for cancer therapy [71].



Figure 13. Distribution of the rate of temperature increase of the agents in vivo.



Figure 15 compares the agent's temperatures when irradiated in vitro and in vivo. This graph shows the differences in heating between the two conditions. This data is essential for understanding how the affects biological environment the agents' performance. It also helps assess whether laboratory results can be expected to be applied to practical applications in the human or animal body. A difference is observed in the agent's temperature when exposed to in vivo radiation, generally resulting in lower temperatures than when the agent is irradiated in an in vitro environment. This change may be due to the laser energy absorbed by the medium (experimental animal) before reaching the agent. However, some agents exhibit higher temperatures when tested in vivo than in vitro. For example, DPP-BT NP exhibited a temperature of 43.6°C at 100 µg/mL during in vitro testing and reached 54°C at 1 mg/mL in 200 µL in vivo [16]. Furthermore, the ISSzyme agent reached a temperature of 54.17°C when irradiated at a wavelength of 808 nm (1 W/cm², 600 s) in vitro and

Conclusion

radiation source.

Components used as therapeutic and contrast agents in photoacoustic imaging-guided PTT have been extensively developed, displaying variations in organic, inorganic, and hybrid organic-inorganic compounds. However, further development of dye-type agents is required. Generally, therapeutic and contrast agent development tends to exhibit absorption in the NIR-I region compared with that in the NIR-II region. Developing agents with NIR-II absorption can offer advantages, including improvement in the contrast of photoacoustic imaging. These agents were then tested for PTT in vitro and in vivo, with 4T1 cells being the most frequently used. When irradiated with specific power and duration, the agents experience a temperature increase of over 40°C, which is sufficient to eliminate cancer cells. No specific relationship was found between the photothermal conversion efficiency and the temperature increase achieved by the agents. Enhancing photothermal conversion not only facilitates improved efficiency in PTT but also elevates photoacoustic signals. Therefore, a deeper exploration of the influence of photothermal conversion efficiency from other perspectives is necessary. There is substantial potential for extensive research involving further exploration in developing agents that act as therapeutic agents and serve as contrast agents to enhance photoacoustic imaging. The development of agents in the NIR-II region holds significant potential because of limited research in this domain.

Future works need to be carried out to develop potential contrast agents that enhance image contrast in photoacoustic imaging and possess therapeutic solid properties. The development of multifunctional agents capable of simultaneously detecting, mapping, and treating diseases will be a primary focus of future research. Additionally, there is potential to improve photothermal conversion efficiency of contrast agents and achieve more uniform heat distribution. Future research should focus on developing new materials with optimized optical and thermal characteristics to enhance the efficiency of photothermal therapy. Although there has been a significant amount of promising preclinical research, broader clinical validation is still needed to ensure the safety and effectiveness of these agents in photoacoustic and photothermal therapy for human treatment.

Further clinical studies will provide a better understanding of the clinical potential of these techniques and accelerate their application in clinical practice. Additionally, challenges remain in translating photoacoustic imaging and photothermal therapy techniques from the laboratory to the clinical practice. Future research should focus on developing methods that can facilitate the integration of these techniques into everyday clinical practice, including the development of more affordable and user-friendly hardware and the creation of standardized clinical protocols. Integrating photoacoustic imaging with other imaging modalities (e.g., MRI, CT) and combining photothermal therapy with other treatment strategies (e.g., chemotherapy, immunotherapy) could provide more comprehensive diagnostic and therapeutic solutions. Therefore, future research is expected to significantly contribute to developing more effective and safe photoacoustic and photothermal therapies.

Table 4. Summar	y of the photo	thermal effects	of the agents in vivo
-----------------	----------------	-----------------	-----------------------

AgentsCellsTherapyToolsAuto(m)Proc Mark AProc Mark AProc (N)BDPAAANANGMPSA11PTPTPLB E35.%31%31%1 </th <th></th>										
PIA-ADMENTERN IT PIA PIA PIA PIA <th< th=""><th>Agents</th><th>Cells</th><th>Therapy</th><th>Tools</th><th>λ_{Rad} (nm)</th><th>Power (W/cm²)</th><th>Time Rad (s)</th><th>ΔT (°C)</th><th>T_{max} (°C)</th><th>Ref</th></th<>	Agents	Cells	Therapy	Tools	λ_{Rad} (nm)	Power (W/cm ²)	Time Rad (s)	ΔT (°C)	T _{max} (°C)	Ref
DPT-PTNPIndIndPTTFUREPSVIR1000001000100100101010101010101010101010101011011 <th< td=""><td>PDA-AuNRs@MSN</td><td>4T1</td><td>PTT</td><td>n/a</td><td>808</td><td>1</td><td>600</td><td>n/a</td><td>61.6</td><td>[15]</td></th<>	PDA-AuNRs@MSN	4T1	PTT	n/a	808	1	600	n/a	61.6	[15]
mHEC-FIEL-AMPKEMEPIC-AMPKVET.VET.VEX.410.55, FLIRE 50.880.810.900<	DPP-BT NP	HeLa	PDT; PTT	FLIR E50;	730	1	600	n/a	54	[16]
ALMUAL MA UK/GC PT n/a 988 1 n/a 988 1 n/a 981 19 51 19 PANLBSKAOT MCF-7 PTT IK termal imaging camera 808; 104 1 800 9.7 4.14 201 PANLBSKAOT MCF-7 (LD) PT IK termal imaging camera 808; 104 1 800 9.7 4.51 26.7 21.1 IBUT-1512 MCF-7 (LD) PT IK termal camera 808 1 300 1.5 5.92 201 IMDS-176_CA HepC2 PTT Rikemal camera 808 1 300 1.5 5.92 206 MIN-172 HepC2 PTT Rikemal camera 603 1.5 300 2.5 5.5 1.25 MIN-172 HepC3 PTT Rikemal camera 608 1 0.00 7.6 6.2 2.01 7.04 COF-560 MTA PTT Rikemal camera 606 1.0.5	mPEG-PEI-AuNRs; mPEG-PEI/CaNPs	MCF-7	UST; PTT	838A-H-O-S; FLIR E5	808	0.8	360; 600	n/a	51	[17]
BITD 302HTPTn/aya <t< td=""><td>Au-MUA-TMA</td><td>U87MG</td><td>PTT</td><td>n/a</td><td>808</td><td>1</td><td>600</td><td>7</td><td>n/a</td><td>[18]</td></t<>	Au-MUA-TMA	U87MG	PTT	n/a	808	1	600	7	n/a	[18]
PANLSRAOTMCPTRRNS	BBTD-1302	4T1	PTT	n/a	980	1	n/a	19	51	[19]
1TMdS, samadnis, 2H-MdS;AS9PTIR hermal langing camera808, 10444025, 31, 265, -721BDF-FR DN-NHeLaPTn/a70n/a601.5121BDF-FNNNHeLaPTn/a70n/a601.5121BDF-FNNNHelp-2PTR heural camera80813001.71.71.7HINCS-FDG-GAHepC2PTTR heural camera80813001.69.72.51.6FAC-SCOMDA-MD21PTTR heural camera80816001.69.89.81.60.01.69.81.60.01.69.81.60.01.69.81.60.01.65.51.21.6 <td>PANI-ES@AOT</td> <td>MCF-7</td> <td>PTT</td> <td>IR camera</td> <td>808: 1064</td> <td>1</td> <td>300</td> <td>9.47</td> <td>46.14</td> <td>[20]</td>	PANI-ES@AOT	MCF-7	PTT	IR camera	808: 1064	1	300	9.47	46.14	[20]
IDT-FR IDT-F IDI-PT IB (DF-FA) 7.0	1T-MoS ₂ nanodots: 2H-MoS ₂	A549	PTT	IR thermal imaging camera	808: 1064	1	480	28.3:12	65.7: -	[21]
IDPA-INF. Hela PTT n/a Name 730 n/a 600 16.9 54 [23] AMDS9815A-NYD Hep-G2 PTT Rikmaging camera 808 1 300 14.7 4.7.7 [25] MNP12 Hep-G2 PTT Rikmaging camera 808 1 300 14.7 4.7.6 [25] PTDPP, PTTPPSNAP MC HC-J PTT Rikmaging camera 808 1 600 2.5 5.9.2 [28] COF-S66 NF HT PDT; PTT n/a 808 1 600 1.6 1.0 [28] COC 5.6 [20] COF-S66 NF HT PTT<	BDT-TH: BDT-NDI	MCF-7	(LED) PTT	IR thermal camera	660	02	600	n/a	53	[22]
ALL DSB SA:N: NYs Hep-2 PTT Flue F100 608 1 300 17.5 7.4 7.4 HMCS: PEG-GA Hep-2 PTT IR imaging camera 808 1 300 1.7 47 [2] HMCS: PEG-GA Hep-2 PTT IR thermal camera 808 1 300 1.7 5.2 [2] FLACS: GO MDA-MB-231 PTT IR thermal camera 635 1.5 300 1.6 5.2 [2] COF-366 NFs 4T1 PTT Rule T100 1064 1.8 0.0 1.7.6 5.2 [3] SQ-NFs 4T1 PTT Rule T400 1.064 1.8 0.0 1.4 5.4 [3] SQ-NFs 4T1 PTT Rule T400 808 1 300 1.9 5.1 [3] SN080 NFs 4T1 PTT Fuke T400 808 1 300 1.4 5.0 1.6 5.0 1.6 5.0 1.7 5.1 5.0 1.6 5.0 1.6 5.0 1.7 1.5	BDP-Fe NPs	HeLa	PTT	n/a	730	n/a	600	16.9	54	[23]
$ \begin{array}{cccccccccccccccccccccccccccccccccccc$	AB-DS@BSA-N2 NYs	Hen-?	PTT	Fluke Ti400	808	1	300	17.5	01	[24]
$ \begin{array}{cccccccccccccccccccccccccccccccccccc$	HMCS-PEG-GA	HepG2	PTT	IR imaging camera	808	1	300	14.7	43.7	[25]
$ \begin{array}{cccccccccccccccccccccccccccccccccccc$	MNPH2	Hop 2	PTT	IR thormal camera	808	1	300	25	50.2	[26]
$\begin{aligned} & \text{PRC-PC} & \text{PTTOPP}, \text{PTOPP}, \text{PTTOPP}, \text{PTOPP}, \text$	FA CS CO	MDA MB 231	PTT	IR thermal camera	808	2	300	2.5 m/a	57.6	[20]
PTIDTP_FTIDTP_SNAPNTS MCP2 PTID	DETORD, DETORD CNIAD NR-	MCE 7	DTT		808	1	500	11/a 20	0	[20]
COP SNYS 411 PTO, PT1 INMERTIAL Gamera 0.5 1.5 5.0 0.7 5.1 2.5 5.2 2.3 5.3 5.0 0.7 5.5 5.2 [33] SQN PN 411 PTT Floke T400 1064 1.8 300 1.7 5.5 [32] GROB PT NPS C6, CL261 PTT Ref Hermal camera 1064 1.8 300 1.9 51 [35] SSAMP 411 PTT IR thermal camera 808 1.0 600 n/a 52 [37] POD-PEG 411 PTT IR thermal camera 808 1.0 600 n/a 60 [38] [34] [37] [37] [37] [37] [36]	COE 2((ND-	MCF-7	F I I DDT. DTT	II/a IB the sum of some sum	606 (25	1	200	20	6U E1	[20]
CA Mank composite 411 PT1 PLORE 1480 868 1 800 1.0 80.2 1.20 molybdenum(V)-based POMs 4T1 PTT Floke T400 1064 1.8 300 2.3 5.5. [32] GCOWET NPs CG.CL261 PT1 R thermal camera 1064 1.6.5 300 n/a 62 [34] SYSTMON NPs 4T1 PTT R tack T400 808 1 600 n/a 51 [35] POMe-NPs 4T1 PTT IK thermal camera 808 1 600 n/a 52 [37] PIG-PTP-GRONS 4T1 PTT IK thermal camera 808 1 300 1.8 6.0 n/a 62.9 [39] POS-PTA NPs 4T1 PTT IR thermal camera 808 1 600 n/a 6.0 1.4 6.0 1.4 6.0 1.4 1.4 PTS-rPANS 4T1 PTT Infared thermography 808	CDA representation	411	PDI; PII	Elsilia Ti490	000	1.5	500	n/a 17.0	51	[29]
SQL NPS 411 PT1 P1/a P1/a <t< td=""><td>CPA nanocomposite</td><td>411</td><td>PII</td><td>Fluke 11480</td><td>808</td><td>1</td><td>600</td><td>17.6</td><td>56.2</td><td>[30]</td></t<>	CPA nanocomposite	411	PII	Fluke 11480	808	1	600	17.6	56.2	[30]
maynetic mody lossed POMs 411 P11 Place 1400 104 1.8 300 2.5 52.5 122 GGD0PT PNs CG CL261 PT IR kardiation Thermal camera 1064 1.50 300 n/a 54 [34] Styme 4T1 PT IR kardiation Thermal camera 808 1 600 n/a 51 [36] BOMs-NPs 4T1 PDT; PTT n/a Rude 808 1 600 n/a 52 [37] The-P-PEG NPs 4T1 PDT; PTT n/a 808 1 800 n/a 60 [37] The-P-PEG NPs 4T1 PTT IR thermal camera 808 1 600 n/a 62 [41] PCS-TPA NPs 4T1 PTT NIR camera 808 1.2 300 n/a 63 640 1/a 631 [45] PCS-TPA NPs 4T1 PTT NIR camera 808 1.4 600 n/a 63	SQ-NPs	411	PII	n/a	760	0.6	900	12	F0 F	[31]
CKCJ@P1 NPs Cb; CL2b1 PT1 IR thermal camera 1064 1; 0.5 300 n/a 54 [34] SSSyme 4T1 PTT IR tadiation Thermal camera 808 1 300 19 51 [36] SV1080 NPs 4T1 PTT IR thermal camera 808 1 600 n/a 52 [37] N16-Pr-PEG NPs 4T1 PTT IR thermal camera 808 1 600 n/a 52.9 [39] N16-Pr-PEG NPs 4T1 PTT IR thermal camera 808 0.8 600 1.8 60 1.4 60 [41] PTA NPs 4T1 PTT IR thermal camera 808 1.2 600 n/a 62 [41] PT, and T AT1 PTT Infrared thermography 808 1.5 600 n/a 63 [43] PGs/ated Sb-835 NPHs 4T1 PTT n/a 630 1.5 600 1.4 63 [44] AudMAs-Cy7 NCs 4T1 PTT Infrared thermography 808 <t< td=""><td>molybdenum(V)-based POMs</td><td>411</td><td>PII</td><td>Fluke 11400</td><td>1064</td><td>1.8</td><td>300</td><td>23</td><td>53.5</td><td>[32]</td></t<>	molybdenum(V)-based POMs	411	PII	Fluke 11400	1064	1.8	300	23	53.5	[32]
Essryme 411 PTI Ikradiaton Inermal camera 808 2 6.00 n/a 6.2 [35] SV1080 NPs 4TI PDT; PTT n/a 755 1 6.00 n/a 51 [35] B <ome-nps< td=""> 4TI PDT; PTT n/a 808 1 6.00 n/a 52 [37] B<ome-nps< td=""> 4TI PDT; PTT n/a 808 1 6.00 n/a 52 [37] J16-PT-PEG NPs 4TI PTT IR thermal camera 808 1 6.00 n/a 60 [41] POS-TPA NPs 4TI PTT IR thermal camera 808 1.5 6.00 n/a 6.2 [40] Mel-PDA-ICG 4TI PTT IR amera 808 1.5 6.00 n/a 6.3 [43] PCS-TPA NPs 4TI PTT n/a 808 1.6 6.00 n/a 6.3 [43] PLG-AtaG 4TI PTT n/a 808 1.6 6.00 n/a 5.2 [40] <</ome-nps<></ome-nps<>	cRGD@PT NPs	C6; GL261	PIT	IR thermal camera	1064	1; 0.5	300	n/a	54	[33]
SY1080 N1/5 411 PTI Fluke T400 808 1 300 19 51 [55] BOMB-NPS 4T1 PTI IR thermal camera 808 1 600 n/a 52 [37] PGP-PEG 4T1 PTT IR thermal camera 808 1 600 n/a 62 [40] PQs-FPA NPS 4T1 PTT IR thermal camera 808 1.2 300 n/a 62 [40] PQs-FPA NPS 4T1 PTT IR thermal camera 808 1.2 500 n/a 62 [41] PT, PT 4T1 PTT IR thermal camera 808 1.5 600 n/a [53] [54] [34] PGS-IPA NPS 4T1 PTT n/a 808 1.6 600 n/a [51] [33] PGS-IPA NPS 4T1 PTT n/a Rameral camera 808 1.6 600 n/a [34] [34] PLACSQD </td <td>ISSzyme</td> <td>411</td> <td>PIT</td> <td>IR radiation Thermal camera</td> <td>808</td> <td>2</td> <td>600</td> <td>n/a</td> <td>62</td> <td>[34]</td>	ISSzyme	411	PIT	IR radiation Thermal camera	808	2	600	n/a	62	[34]
B-OMe.NPs 4F1 PDT, PTT n/a 735 1 600 n/a 51 [56] QDas-PEG 4T1 PDT; PTT n/a 808 1 600 n/a 52 [37] Y16-Pr-PEG NPs 4T1 PTT IR thermal camera 808 1 600 n/a 60 [39] QPs-TPA NPs 4T1 PTT IR thermal camera 808 1.2 300 n/a 60 [41] PT, m/T 4T1 PTT NIR camera 808 1.5 600 n/a [42] M1 NPs 4T1 PTT n/a 808 1.5 600 n/a [43] J-Au-CSQ 4T1 PTT n/a 808 1.5 600 n/a [42] J-Au-CSQ 4T1 PTT IR thermal camera 808 1.0 600 n/a [42] J-Au-CSQ 4T1 PTT IR thermal camera 808 1.0 500 n/a </td <td>SY1080 NPs</td> <td>4T1</td> <td>PTT</td> <td>Fluke Ti400</td> <td>808</td> <td>1</td> <td>300</td> <td>19</td> <td>51</td> <td>[35]</td>	SY1080 NPs	4T1	PTT	Fluke Ti400	808	1	300	19	51	[35]
BQDs-FIG 4F1 PT IR thermal camera 808 1 600 n/a 52 [37] VIG-P-FEG NPs 4T1 PTT IR thermal camera 808 1 600 n/a 602 [38] AuNSPHs 4T1 PTT IR thermal camera 808 1.2 300 n/a 62 [40] Mel-pDA-ICG 4T1 PTT IR thermal camera 808 1.2 600 n/a 62 [43] PEGylated Sb-835 NPHs 4T1 PTT Infared thermography 808 1 600 n/a 63.3 [43] PEGylated Sb-835 NPHs 4T1 PTT n/a 808 1.5 600 n/a 63.3 [45] JAu-CSQD 4T1 PTT n/a 808 1.6 600 n/a 63.3 [44] AuuMbA_CY7 NCs 4T1 PTT Rubernal camera 808 1.3 61.4 [48] VAuCSQ KOV-3 PTT <td< td=""><td>B-OMe-NPs</td><td>411</td><td>PDT; PTT</td><td>n/a</td><td>735</td><td>1</td><td>600</td><td>n/a</td><td>51</td><td>[36]</td></td<>	B-OMe-NPs	411	PDT; PTT	n/a	735	1	600	n/a	51	[36]
Y16-Pr-PEGNPs 4I1 PDT, PIT n/a 808 1 600 n/a 60 [38] AuNSPHs 4I1 PIT R thermal camera 808 0.8 0.00 2.8 62 [40] Mel-pDA-ICG 4I1 PIT NIR camera 808 1.2 300 n/a 60 [41] PIT; mPT 4I1 PIT SolidSpec-3700 808 1.6 600 n/a 543 PEGylatedSb-83SNPHs 4I1 PIT n/area 808 1 600 n/a 61.3 [41] AuMBA _{2x} -Cy7NCs 4T1 PIT n/a 808 1 600 n/a 61.3 [44] AuABA _{2x} -Cy7NCs 4T1 PIT n/a Rand 808 1 600 n/a 61.3 [44] AuABA _{2x} -Cy7NCs 4T1 PIT IR thermal camera 808 1.6 600 1.3 60.4 1.4 61.3 61.3 61.4 61.3 61.3 61.4 61.3 61.3 61.4 61.3 61.3 61.4 </td <td>BQDs-PEG</td> <td>4T1</td> <td>PTT</td> <td>IR thermal camera</td> <td>808</td> <td>1</td> <td>600</td> <td>n/a</td> <td>52</td> <td>[37]</td>	BQDs-PEG	4T1	PTT	IR thermal camera	808	1	600	n/a	52	[37]
AuNSPHS 4T1 PTT IR thermal camera 808 1 300 18.9 52.9 [39] PQS-TPA NPS 4T1 PTT IR thermal camera 808 1.2 300 n/a 60 [41] P1: mT 4T1 PTT SolidSpec-3700 808 1.2 300 n/a 6.0 [41] P1: mT 4T1 PTT Infrared thermography 808 1.5 600 n/a 5.5 [43] P40 AudMBAs-Cy NCs 4T1 PTT n/ared thermography 808 1.6 600 n/a 5.2 [46] AuuSAs-Cy NCs 4T1 PTT n/ared thermography 808 1.6 600 n/a 5.2 [46] AuuMBAs-Cy NCs 4T1 PTT Rthermal camera 808 1.6 0.0 n/a 5.2 [46] VauSO-anarods 4T1 CDT; PTT FLR SC300 1061 0.25 0.0 n/a 5.1 [50] VB-So anarods 4T1 CDT; PTT FLR AS20 808 1.8 300 n/a </td <td>Y16-Pr-PEG NPs</td> <td>4T1</td> <td>PDT; PTT</td> <td>n/a</td> <td>808</td> <td>1</td> <td>600</td> <td>n/a</td> <td>60</td> <td>[38]</td>	Y16-Pr-PEG NPs	4T1	PDT; PTT	n/a	808	1	600	n/a	60	[38]
PQs-TPA NPs 4T1 PTT IR kramal camera 808 0.8 600 28 62 [40] PRI-PDA-ICG 4T1 PTT NIR camera 808 1.2 600 n/a 60 [41] PT; mPT 4T1 PTT Infrared thermography 808 1 600 n/a 55 [43] PEOylated 5b-855 NPIs 4T1 PTT n/a 808 1 600 n/a 61.3 [45] J-Au-CS-QD 4T1 PTT n/a 808 1.6 600 n/a 5.2 [46] MAMBA ₃₅ -Cy7NCs 4T1 PTT IR thermal camera 808 1.6 600 n/a 5.2 [47] MNFs 4T1 PTT IR thermal imaging camera 808 1.8 3.00 n/a 4.2 [49] PD-FA nanoparticles 4T1 CDr, PTT Fukersal imaging camera 845 1.8 3.00 n/a 5.1 [56] [50] [50]	AuNSPHs	4T1	PTT	IR thermal camera	808	1	300	18.9	52.9	[39]
Mel-pDA-ICG 4T1 PTT NIR camera 808 1.2 300 n/a 600 [41] PT; mPT 4T1 PTT SolidSpec-3700 808 1.5 600 n/a 55 [43] PIGylated 5b-835 NPHs 4T1 PTT Infrared thermography 808 1 600 n/a 55 [43] PIGylated 5b-835 NPHs 4T1 PTT n/a mode 808 1 600 n/a 61.3 [45] Au_MBA_S-Cy7NCs 4T1 PTT IR thermal camera 808 1.6 600 n/a 61.3 [45] JAu-CS-QD 4T1 PTT Fluke Ti400 808 1.8 30.0 n/a 62.4 [49] PD-FA nanoparticles HeLa PTT Fluke Ti400 808 0.8 300 n/a 65 [51] Fe-ZD5 PPT-Sc IPP-Sc IPP-Sc IPP-Sc NPs A549 PTT Ithermal imager 808 0.82 360 n/a 55 [54] PDP-Sc IPP-Sc IPP-Sc IPP-Sc IPP Sc MP A549 PTT Ithermal image	PQs-TPA NPs	4T1	PTT	IR thermal camera	808	0.8	600	28	62	[40]
PT; mT4T1PTSolidspec-37008081.5600n/am/a[42]M1 NPs4T1PTTInfrared hermography808154058n/a[44]M1 MPs4T1PTTn/a808154058n/a[44]AuuMBAas-Cy7 NCs4T1PTTn/a8081600n/a61.3[45]J-Au-CS-QD4T1PTTR thermal cmera8080.5600n/a61.3[46]GCO-AuNPSKOV-3PTTFLIR SC30010610.25300n/a61.4[47]MNFs4T1PTTFLIR SC3008081300n/a[51][51]PD-FA nanoparticlesHeLaPTTFLIR A3008081300n/a[52][51]PD-FA nanoparticlesHeLaPTTIthermal imaging camera8451300n/a[52][51]Fe-ZDSEMT-6PTTIthermal imagera8080.8300n/a[52][51]PD-FS: DPP-SE: NP-SEN PSA549PTTIthermal imagera8080.8300n/a[54][54]PD-FS: DPP-SE: NP-SEN PSA549PTTIthermal imagera8081.26001.83.6[54]PD-FS: DPP-SE: NPS-SEN PSA549PTTIthermal imagera8081.26001.83.6[54]PDF-SS: DPP-SE: NPS-SEN PSA549PTTIthermal	Mel-pDA-ICG	4T1	PTT	NIR camera	808	1.2	300	n/a	60	[41]
M1NsM1M1M1M1m/a8081600n/a55[43]PEGylated Sb-835 NPHs411PTn/a8081600n/a61.3[41]AuuMBAa-Cy7 NCs411PTn/a80816001.361.3[41]J-Au-CS-QD411PTR thermal camera8081.66000.362.8[41]GCA-UNPSKOV-3PTFILR SC30010610.25300n/a61.4[41]WIS-Go nanorods411PTTFILR A3008081.0300n/a62.1[51]PD-FA nanoparticlesHcl.aPTTFILR A3008080.830.0n/a65[51]PD-FSC; DPP-SS; DPP-SSe NPsA549PTTFluke Ti4808080.830.0n/a65[51]Fe-ZDSEMT-6PTTR thermal imager8080.830.0n/a65[51]PD-Sc; DPP-SF; DPP-SF; DPP-SFA549PTTR thermal imager8080.80.0n/a55[51]Fe-ZDSA549PTTR thermal imager8080.80.0n/a65[51]DPP-Sc; DPP-SF; DPP-SF; DPP-SFA549PTTR thermal imager8080.80.0n/a56.5[51]IDT-BTZMDA-TEGA549PTTR tamera8081.06001.450.2[51]DPP-Sc; DPP-SF; DPP-SF; DPP-SFM1<	PT; mPT	4T1	PTT	SolidSpec-3700	808	1.5	600	n/a	n/a	[42]
PECylated Sb-835 NPHs fTI PIT n/a Na	M1 NPs	4T1	PTT	Infrared thermography	808	1	600	n/a	55	[43]
Aua,MBab-Cy7 NCs 4T1 PTT I/a Na 808 1 600 n/a 61.3 [45] J-Au-CSQD 4T1 PTT IR thermal camera 808 0.5 600 .04 52.8 [46] J-Au-CSQD SKOV-3 PTT IR thermal camera 808 1 600 .04 52.8 [47] MNFs 4T1 DTF FLIK SC300 808 1 600 .08 51.4 [48] WisOp nanorods 4T1 DTF FLIK A300 808 1 .000 .04 42 [49] PD-FA nanoparticles ICL ICL Rthermal imaging camera 845 1 .000 .0	PEGylated Sb-835 NPHs	4T1	PTT	n/a	808	1	540	58	n/a	[44]
J-Au-CS-QD 4T1 PT IR thermal camera 808 0.5 600 30.4 52.8 [4] rGO-AUNP SKOV-3 PTT FLIR SC300 1061 0.55 300 n/a n/a 1/a MNFs 4T1 CDT; PTT FLIR A300 808 1 300 n/a 42 [49] PD-FA nanoparticles HeLa CDT; PTT IR thermal inaging camera (CEM, at an interval of 60 s) 808 1.8 300 n/a 65 [51] PD-FS, DPP-SS, DPP-SS NPS A549 PTT IR thermal imaging camera (CEM, at an interval of 60 s) 808 0.82 360 n/a 65 [51] Fe-ZDS EMT-6 PTT IR thermal imager 808 0.88 0.00 n/a 54 [54] PDM tumor CDT; PTT IR thermal imager 808 0.75 600; 300 n/a 55 [54] PDP-SS, DPP-SF; DPP-SF A549 PTT Rice Ti480 808 1.2 60; 300 1.4 56 56.7 [55] TMB-F4TCNQ; TMB-TCNQ 4	Au ₄₄ MBA ₂₆ -Cy7 NCs	4T1	PTT	n/a	808	1	600	n/a	61.3	[45]
rGO-AuNP SKOV-3 PIT FLIR SC300 1061 0.25 300 n/a n/a [4] MNFs 4T1 PIT FLIR A100 808 1 600 18.3 51.4 [4] WisQue nanorods HELa CDT; PTT FLIR A300 808 1.0 300 n/a 42 [4] PD-FA nanoparticles HELa CDT; PTT Rthermal imaging camera (CEM, at an interval of 60 s) 885 0.8 300 n/a 65 [5] DP-So; DPP-Ss; DPP-Se NFS A549 PTT Iktermal imager 808 0.82 300 n/a 612 [5] POM tumor CDT; PTT Iktermal imager 1060 0.82 600 n/a 7.4 [5] DPP-Ss; DPP-S; DPP-SeF NP A549 PTT Iktermal imager 1060 0.82 300 n/a 5.6 [5] IDCL NPs Af10 PTT Iktermal imaging camera 808 1.2 600.0 1.4 30.2 5.6 [5] IDCL NPs 411 PTT Ikter Ti400	J-Au-CS-QD	4T1	PTT	IR thermal camera	808	0.5	600	30.4	52.8	[46]
MNFs4T1PTTFluke Ti400808160018.351.4[4]WisOs nanorods4T1CDT; PTTFLIR A3008081300n/a42[4]PD-FA nanoparticlesHEAPTTIR thermal imaging camea845300n/a650[5]DPF-SC; DPP-SS; DPP-SS NPSA549PTTFluke Ti4808080.82300n/a65[5]Fe-ZDSMT-6PTTIR thermal imager8080.82300n/a52[5]PD-S5; DPP-SF; DPP-SF NPA549PTTIR thermal imager8080.82300n/a55[5]PD-S5; DPP-OF; DPP-SF; DPP-SF NPA549PTTIR thermal imager8080.72300n/a55[5]DPT-S7A549PTTIR thermal imager8081.2600; 300n/a55[5]IDCL NPA549PTTIR thermal imager8081.2600; 3001.656[5]DCL NPS4T1PTTIR thermal camera8081.2600; 3001.656[5]IDCL NPS4T1PTTIL thermal camera8081.160024.354.5[6]IDCL NPS4T1PTTIL thermal imaging camera8081.1300n/a58.5[6]AuGu_tNPS ChAPS4T1PTTIL thermal imaging camera8081.1300n/a50.5[6]AuGu_tNPS ChAPS <td>rGO-AuNP</td> <td>SKOV-3</td> <td>PTT</td> <td>FLIR SC300</td> <td>1061</td> <td>0.25</td> <td>300</td> <td>n/a</td> <td>n/a</td> <td>[47]</td>	rGO-AuNP	SKOV-3	PTT	FLIR SC300	1061	0.25	300	n/a	n/a	[47]
WisOse nanorods 4T1 CDT; PTT FLIR A300 808 1 300 n/a 42 [4] PD-FA nanoparticles HeLa PTT IR thermal imaging camera (CEM, at an interval of 60 s) 845 1 300 n/a 65 [5] DPP-SO; DPP-SS; DPP-SSe NPs A549 PTT Fluke Ti480 808 0.8 300 n/a 65 [5] POM tumor CDT; PTT IR thermal imagers 808 0.82 360 n/a 65 [5] POM tumor CDT; PTT IR thermal imagers 808 0.82 360 n/a 65 [5] POM tumor CDT; PTT IR thermal imagers 808 1.2 600, 0 n/a 55 [5] PDF-Sc; DPP-SF; DPP-SF A549 PTT Fluke Ti480 808 1.2 600, 0 1.4 53 56 DTGC NPs 4T1 PTT FLIR A50 630 1.4 600 4.5 6.2 [5] <td>MNFs</td> <td>4T1</td> <td>PTT</td> <td>Fluke Ti400</td> <td>808</td> <td>1</td> <td>600</td> <td>18.3</td> <td>51.4</td> <td>[48]</td>	MNFs	4T1	PTT	Fluke Ti400	808	1	600	18.3	51.4	[48]
PD-FA nanoparticles HeLa PTT IR thermal imaging camera (CEM, at a interval of 60 s) 845 1 300 n/a 65 [50] DPP-SQ; DPP-SS; DPP-SSe NPs A549 PTT Fluke Ti480 808 0.82 300 n/a 65 [51] Fe-ZDS EMT-6 PTT IR thermal imager 808 0.82 300 n/a 65 [51] POM tumor CDT; PTT IR thermal imager 1060 0.8 600 n/a 65 [54] DPP-SS; DPP-OF; DPP-SF; DPP-SeFN A549 PTT Fluke Ti480 808 0.75 300 n/a 55 [54] MPDA-PEG A549 PTT IR imager 600 1 000 18.8 3.6 [56] IDCIC NPS 4T1 PTT Fluke Ti400 808 1 000 1.8 3.6 [56] IDCIC NPS 4T1 PTT Thermal camera 808 1.1 50 0.0 n/a 5.4 6.3 </td <td>W₁₈O₄₉ nanorods</td> <td>4T1</td> <td>CDT; PTT</td> <td>FLIR A300</td> <td>808</td> <td>1</td> <td>300</td> <td>n/a</td> <td>42</td> <td>[49]</td>	W ₁₈ O ₄₉ nanorods	4T1	CDT; PTT	FLIR A300	808	1	300	n/a	42	[49]
DPP-SQ: DPP-SQ: DPP-SQ: PNPSA549PTFurce T4808080.8300n/a65[51]Fe-ZDSEMT-60PTTIR thermal imager8080.8236045n/a[52]POMtumorCDT; PTIR thermal images10600.8600n/a[53]DPP-SQ: DPP-SF; DPP-SF PNA549PTFluke Ti4808080.75300n/a55[54]MPDA-PEGA549PTIR camera8081,260; 30026.656.7[55]TMB-F4TCNQ; TMB-TCNQ4T1PTIR imager1060130018.853.6[56]DCIC NPsGBP-FANBHepG2; PC3PTFLIR Ax560160024.5060.2[57]Au@Cu ₂ -xSe NPs4T1PTThermal camera808160024.5060.2[57]Au@Cu ₂ -xSe NPs4T1PTThermal camera8081.860085.5[59]Au@Cu ₂ -xSe NPs4T1PTThermal imaging camera8081.1300n/a50.1[61]MnCo@CNPpMCF-7PTThermal imaging camera8081.830.0n/a53.5[64]GNS; Silica@GNSMDA-MB-21PTInker Ti5S5660; 8080.76; 2300n/a53.5[64]GNR@MNPs-Gd-64 CuHep-2PTFluke Ti5D606; 8081.5300n/a53.5[64]GNR@MN	PD-FA nanoparticles	HeLa	PTT	IR thermal imaging camera	845	1	300	n/a	65	[50]
DPP-SS; DPP-SS; DPP-SS NPS A549 P11 Fluke 1480 808 0.8 300 n/a 65 [51] Fe-ZDS EMT-6 PTT Ik thermal imager 808 0.82 360 45 n/a [52] POM tumor CDT; PTT IR thermal images 1060 0.8 600 n/a n/a [53] DPP-SS; DPP-SF; DPP-SF; DPP-SF N A549 PTT Fluke Ti480 808 0.75 300 n/a 55 [54] MPDA-PEG A549 PTT IR camera 808 1;2 600;300 26.6 56.7 [55] TMB-F4TCNQ; TMB-TCNQ 4T1 PTT IR imager 1060 1 300 18.8 53.6 [56] DCIC NPs 4T1 PTT Fluke Ti400 808 1 600 24.63 54.63 [58] IDT-BTZTD Pdots 4T1 PTT Thermal camera 808 0.8 0.8 0.60 58.5 [59] Au@CQ_2xSe NPs 4T1 PTT Thermal imaging camera 808 1 30				(CEM, at an interval of 60 s)				,	-	[=-1]
Fe-ZDS EM1-6 PTT IR thermal imager 808 0.82 360 45 n/a [52] POM tumor CDT; PTT IR thermal images 1060 0.8 600 n/a n/a [53] DPP-SS; DPP-OF; DPP-SF; DPP-SF A549 PTT IR thermal imager 808 0.75 300 n/a 55 [54] MPDA-PEG A549 PTT IR camera 808 1.2 600; 300 26.6 56.7 [55] IDCIC NPs 4T1 PTT IR mager 1060 1 300 n/a 53.6 [56] IDCIC NPs 4T1 PTT FILR Ax5 630 1 600 24.63 54.63 [58] IDT-BTZTD Pdots 4T1 PTT Thermal camera 808 0.8 600 58.5 [59] Au@Cu ₂ -xSe NPs 4T1 PTT Thermal camera 808 1 300 n/a 50.1 [60] FA-DOX-ICG-PFP@Lip Y79 PTT Thermal camera 808 1 300 n/a 53.5<	DPP-SO; DPP-SS; DPP-SSe NPs	A549	PIT	Fluke 11480	808	0.8	300	n/a	65	[51]
POM tumor CDT; PTT IR thermal images 1060 0.8 600 n/a n/a [53] DPP-SS; DPP-SF; DPP-SeF NP A549 PTT Fluke Ti480 808 0.75 300 n/a 55 [54] MPDA-PEG A549 PTT IR camera 808 1; 2 600; 300 26.6 56.7 [55] TMB-F4TCNQ; TMB-TCNQ 4T1 PTT IR imager 1060 1 300 18.8 53.6 [66] DCIC NPs 4T1 PTT Fluke Ti400 808 1 600 24.63 54.63 [58] IDT-BTZTD Pdots 4T1 PTT Thermal camera 808 0.8 600 - 58.5 [59] Au@Cu ₂ -xSe NPs 4T1 PTT Thermal camera 808 1 300 n/a 50.1 [61] MnCo@CNPpy Y79 PTT Thermal imaging camera 808 1 300 n/a 53.5 [63] ZPC NPs MDA-MB-231 PTT n/a 1 300 n/a 53.5	Fe-ZDS	EMT-6	PTT	IR thermal imager	808	0.82	360	45	n/a	[52]
DPP-SS; DPP-OF; DPP-SF; DPP-SF; DPP-SF NP A549 PTT Fluke Ti480 808 0.75 300 n/a 55 [54] MPDA-PEG A549 PTT IR camera 808 1;2 600; 300 26.6 56.7 [55] TMB-F4TCNQ; TMB-TCNQ 4T1 PTT IR imager 1060 1 300 18.8 53.6 [56] IDCIC NPs 4T1 PTT Fluke Ti400 808 1 600 24.63 64.63 [57] GBP-FANB HepG2; PC3 PTT FLIR Ax5 630 1 600 24.63 54.63 [59] Au@Cu ₂ -xSe NPs 4T1 PTT Thermal camera 808 0.8 600	POM	tumor	CDT; PTT	IR thermal images	1060	0.8	600	n/a	n/a	[53]
MPDA-PEG A549 PIT IR camera 808 1; 2 600; 300 26.6 56.7 [55] TMB-F4TCNQ; TMB-TCNQ 4T1 PIT IR imager 1060 1 300 18.8 53.6 [56] IDCIC NPs 4T1 PIT Fluke Ti400 808 1 600 24.5 60.2 [57] GBP-FANB HepG2; PC3 PIT FLIR Ax5 630 1 600 24.63 54.63 [58] IDT-BTZTD Pdots 4T1 PIT Thermal camera 808 0.8 600 58.5 [59] Au@Cu ₂ -xSe NPs 4T1 PIT Thermal camera 808 1 300 n/a 50.1 [61] MnCO@CNPpy MCF-7 PIT Temperature transmitter 808 1 300 n/a n/a [62] GNS; Silica@GNS MDA-MB-231 PIT n/a 808 1 900 n/a 53.5 [64] GNR@MNPs-Gd-64 Cu Hep-2 PIT Fluke Ti400 808 1.5 300 n/a 58.5	DPP-SS; DPP-OF; DPP-SF; DPP-SeF NP	A549	PTT	Fluke Ti480	808	0.75	300	n/a	55	[54]
TMB-F4TCNQ: TMB-TCNQ 4T1 PTT IR imager 1060 1 300 18.8 53.6 [56] IDCIC NPs 4T1 PTT Fluke Ti400 808 1 600 24.5 60.2 [57] GBP-FANB HepG2; PC3 PTT FLIR Ax5 630 1 600 24.63 54.63 [58] IDT-BTZTD Pdots 4T1 PTT Thermal camera 808 0.8 600 58.5 [59] Au@Cu ₂ -xSe NPs 4T1 PTT Thermal camera 808 1.1 300 n/a 50.1 [61] MnCO@CNPpy MCF-7 PTT Temperature transmitter 808 1 300 n/a n/a [62] GNS; Silica@GNS MDA-MB-231 PTT n/a 808 1 900 n/a n/a [63] ZPC NPs 4T1 PDT; PTT Fluke Ti555 660; 808 0.76; 2 300 n/a 53.5 [64] GNR@MNPs-Gd-64 Cu Hep-2 PTT Fluke Ti400 808 1 600 n/a 58.	MPDA-PEG	A549	PTT	IR camera	808	1; 2	600; 300	26.6	56.7	[55]
IDCIC NPs 4T1 PTT Fluke Ti400 808 1 600 24.5 60.2 [57] GBP-FANB HepG2; PC3 PTT FLIR Ax5 630 1 600 24.63 54.63 [58] IDT-BTzTD Pdots 4T1 PTT Thermal camera 808 0.8 600 58.5 [59] Au@Cu ₂ -xSe NPs 4T1 PTT Thermal camera 808 1,1.5 300 31.5 48.9 [60] FA-DOX-ICG-PFP@Lip Y79 PTT Thermal imaging camera 808 1 300 n/a 50.1 [61] MnCO@CNPpy MCF-7 PTT Temperature transmitter 808 1 900 n/a n/a [62] GSN Silica@GNS MDA-MB-231 PTT r/a Temperature transmitter 808 1 900 n/a n/a [63] ZPC NPs 4T1 PDT; PTT Fluke Ti555 660; 808 0.76; 2 300 n/a 53.5 [64] GNR@MNPs-Gd-64 Cu Hep-2 PTT Fluke Ti400 808 1.5	TMB-F4TCNQ; TMB-TCNQ	4T1	PTT	IR imager	1060	1	300	18.8	53.6	[56]
GBP-FANB HepG2; PC3 PIT FLIR Ax5 630 1 600 24.63 54.63 [58] IDT-BTZTD Pdots 4T1 PTT Thermal camera 808 0.8 600 58.5 [59] Au@Cu ₂ -xSe NPs 4T1 PTT Thermal camera 808 1; 1.5 300 31.5 48.9 [60] FA-DOX-ICG-PFP@Lip Y79 PTT Thermal imaging camera 808 1 300 n/a 50.1 [61] MnCO@CNPpy MCF-7 PTT Temperature transmitter 808 1 300 n/a n/a [62] GNS; Silica@CNS MDA-MB-231 PTT n/a RS 1.8 900 n/a n/a [63] ZPC NPs 4T1 PDT; PTT Fluke Ti555 660; 808 0.76; 2 300 n/a 53.5 [64] GNR@MNPs-Gd-64 Cu Hep-2 PTT Fluke Ti400 808 1.5 300 n/a 58.5 [65] PFTQ-PEG-Gd NPs 4T1 PTT FLIR E50 808 1 600 n/a<	IDCIC NPs	4T1	PTT	Fluke Ti400	808	1	600	24.5	60.2	[57]
IDT-BTzTD Pdots 4T1 PTT Thermal camera 808 0.8 600 58.5 [59] Au@Cu2-xSe NPs 4T1 PTT Thermal camera 808 1; 1.5 300 31.5 48.9 [60] FA-DOX-ICG-PFP@Lip Y79 PTT Thermal imaging camera 808 1 300 n/a 50.1 [61] MnCO@CNPpy MCF-7 PTT Temperature transmitter 808 1 300 n/a n/a [62] GNS; Silica@CNS MDA-MB-231 PTT n/a 808 1 900 n/a n/a [63] ZPC NPs 4T1 PDT; PTT Fluke Ti555 660; 808 0.76; 2 300 n/a 53.5 [64] GNR@MNPs-Gd-64 Cu Hep-2 PTT Fluke Ti400 808 1.5 300 n/a 58.5 [65] PFTQ-PEG-Gd NPs 4T1 PTT FLIR E50 808 1 600 n/a 58.5 [66] Ag250M-P-RGD HeLa PT IR thermal imaging 808 2 300 n/a	GBP-FANB	HepG2; PC3	PTT	FLIR Ax5	630	1	600	24.63	54.63	[58]
Au@Cu ₂ -xSe NPs 4T1 PTT Thermal camera 808 1; 1.5 300 31.5 48.9 [60] FA-DOX-ICG-PFP@Lip Y79 PTT Thermal imaging camera 808 1 300 n/a 50.1 [61] MnCO@CNPpy MCF-7 PTT Temperature transmitter 808 1 300 n/a n/a [62] GNS; Silica@GNS MDA-MB-231 PTT n/a 808 1 900 n/a n/a [63] ZPC NPs 4T1 PDT; PTT Fluke TiS55 660; 808 0.76; 2 300 n/a 53.5 [64] GNR@MNPs-Gd-64 Cu Hep-2 PTT Fluke Ti400 808 1.5 300 n/a 58.5 [65] PFTQ-PEG-Gd NPs 4T1 PTT FLIR E50 808 1 600 n/a 58.5 [66] Ag_2S@M-P.RGD HeLa PTT IR thermal imaging 808 2 360 n/a 53.3 [67] Au@PB@Cu_20@BPQDs/PAHNCs HeLa PT PTT monitoring system (MG33) 50	IDT-BTzTD Pdots	4T1	PTT	Thermal camera	808	0.8	600		58.5	[59]
FA-DOX-ICG-PFP@Lip Y79 PTT Thermal imaging camera 808 1 300 n/a 50.1 [61] MnCO@CNPpy MCF-7 PTT Temperature transmitter 808 1 300 n/a n/a [62] GNS; Silica@GNS MDA-MB-231 PTT n/a 808 1 900 n/a n/a [63] ZPC NPs 4T1 PDT; PTT Fluke TiS55 660; 808 0.76; 2 300 n/a 53.5 [64] GNR@MNPs-Gd-64 Cu Hep-2 PTT Fluke Ti400 808 1.5 300 n/a 58.5 [65] PFQ-PEG-Gd NPs 4T1 PTT FLIR E50 808 1 600 n/a 58.5 [66] Ag2S@M-P.RGD HeLa PTT IR thermal imaging system (MG33) 608 2 360 n/a 53 [67] Au@PB@Cu_20@BPQDs/PAH NCs HeLa PT PTT monitoring system (MG33) 600 1.5 300 1.6.9 n/a [68]	Au@Cu ₂ -xSe NPs	4T1	PTT	Thermal camera	808	1; 1.5	300	31.5	48.9	[60]
MnCO@CNPpy MCF-7 PTT Temperature transmitter 808 1 300 n/a n/a [62] GNS; Silica@GNS MDA-MB-231 PTT n/a 808 1 900 n/a n/a [63] ZPC NPs 4T1 PDT; PTT Fluke TiS55 660; 808 0.76; 2 300 n/a 53.5 [64] GNR@MNPs-Gd-64 Cu Hep-2 PTT Fluke Ti400 808 1.5 300 n/a 58.5 [65] PFQ-PEG-Gd NPs 4T1 PTT FLIR E50 808 1 600 n/a 58.5 [66] Ag_2S0M-P.RGD HeLa PTT IR thermal imaging 808 2 360 n/a 53 [67] Au@PB@Cu_2O@BPQDs/PAHNCs HeLa PT PTT monitoring system (MG3) 600 1.5 300 n/a 53 [67]	FA-DOX-ICG-PFP@Lip	Y79	PTT	Thermal imaging camera	808	1	300	n/a	50.1	[61]
GNS; Silica@GNS MDA-MB-231 PTT n/a 808 1 900 n/a n/a [63] ZPC NPs 4T1 PDT; PTT Fluke TiS55 660; 808 0.76; 2 300 n/a 53.5 [64] GNR@MNPs-Gd-64 Cu Hep-2 PTT Fluke Ti400 808 1.5 300 n/a 58.5 [65] PFTQ-PEG-Gd NPs 4T1 PTT FLIR E50 808 1 600 n/a 58.5 [66] Ag2S@M-P.RGD HeLa PTT IR thermal imaging 808 2 360 n/a 53 [67] Au@PB@Cu_2O@BPQDs/PAH NCs HeLa PTT PTT monitoring system (MG33) 650 1.5 300 16.9 n/a [68]	MnCO@CNPpy	MCF-7	PTT	Temperature transmitter	808	1	300	n/a	n/a	[62]
ZPC NPs 4T1 PDT; PTT Fluke TiS55 660; 808 0.76; 2 300 n/a 53.5 [64] GNR@MNPs-Gd-64 Cu Hep-2 PTT Fluke Ti400 808 1.5 300 n/a 58.5 [65] PFTQ-PEG-Gd NPs 4T1 PTT FLIR E50 808 1 600 n/a 58.5 [66] Ag_S@M-P-RGD HeLa PTT IR thermal imaging 808 2 360 n/a 53 [67] Au@PB@Cu_20@BPQDs/PAH NCs HeLa PTT PTT monitoring system (MG33) 650 1.5 300 n/a 53 [67]	GNS; Silica@GNS	MDA-MB-231	PTT	n/a	808	1	900	n/a	n/a	[63]
GNR@MNPs-Gd-64 Cu Hep-2 PTT Fluke Ti400 808 1.5 300 n/a 58.5 [65] PFTQ-PEG-Gd NPs 4T1 PTT FLIR E50 808 1 600 n/a 58.5 [66] Ag:S@M-P-RGD HeLa PTT IR thermal imaging 808 2 360 n/a 53 [67] Au@PB@Cu_2O@BPQDs/PAH NCs HeLa PTT PTT monitoring system (MG33) 650 1.5 300 16.9 n/a [68]	ZPC NPs	4T1	PDT; PTT	Fluke TiS55	660; 808	0.76; 2	300	n/a	53.5	[64]
PFTQ-PEG-Gd NPs 4T1 PTT FLIR E50 808 1 600 n/a 58.5 [66] Ag:S@M-P-RGD HeLa PTT IR thermal imaging 808 2 360 n/a 53 [67] Au@PB@Cu_2O@BPQDs/PAH NCs HeLa PTT PTT monitoring system (MG33) 650 1.5 300 16.9 n/a [68]	GNR@MNPs-Gd-64 Cu	Hep-2	PTT	Fluke Ti400	808	1.5	300	n/a	58.5	[65]
AgS@M-P-RGD HeLa PTT IR thermal imaging 808 2 360 n/a 53 [67] Au@PB@Cu ₂ O@BPQDs/PAH NCs HeLa PTT PTT monitoring system (MG33) 650 1.5 300 16.9 n/a [68]	PFTO-PEG-Gd NPs	4T1	PTT	FLIR E50	808	1	600	n/a	58.5	[66]
Au@PB@Cu ₂ O@BPQDs/PAH NCs He La PTT PTT monitoring system (MG33) 650 1.5 300 16.9 n/a [68]	Ag ₂ S@M-P-RGD	HeLa	PTT	IR thermal imaging	808	2	360	n/a	53	[67]
$\Sigma_{\rm c}$, Σ_{\rm	Au@PB@Cu ₂ O@BPODs/PAH NCs	HeLa	PTT	PTT monitoring system (MG33)	650	1.5	300	16.9	n/a	[68]
	PS: n/a signifies that the information pro-	wided is not writte	n or not detail	ad					·,	[~~]



Acknowledgements

We would like to thank Indonesia Endowment Funds for Education (LPDP) and Center for Higher Education Funding (BPPT) for supporting this work.

Author contributions

Donni Kis Apriyanto: data collection, formal analysis, and writing-original draft preparation. Rini Widyaningrum: conceptualization, result validation, and review. Mitrayana: supervision and reviewing. Andreas Setiawan: supervision and reviewing.

Data availability

The evidence data generated and/or analyzed during the current study are available from the corresponding author on reasonable request.

Competing Interests

The authors have declared that no competing interest exists.

References

- 1. Beard P. Biomedical photoacoustic imaging. Interface Focus [Internet]. 2011; 1: 602–31. Available at:
- https://royalsocietypublishing.org/doi/10.1098/rsfs.2011.0028
 Moore C, Jokerst J V. Strategies for Image-Guided Therapy, Surgery, and Drug Delivery Using Photoacoustic Imaging. Theranostics [Internet]. 2019; 9: 1550–71. Available at: http://www.thno.org/v09p1550.htm
- Das D, Sharma A, Rajendran P, Pramanik M. Another decade of photoacoustic imaging. Phys Med Biol [Internet]. 2021; 66: 05TR01. Available at: https://iopscience.iop.org/article/10.1088/1361-6560/abd669

- Das D, Pramanik M. Combined ultrasound and photoacoustic imaging of blood clot during microbubble-assisted sonothrombolysis. J Biomed Opt [Internet]. 2019; 24: 1. Available at: https://www.spiedigitallibrary.org/journals/journal-of-biomedical-optics/v olume-24/issue-12/121902/Combined-ultrasound-and-photoacoustic-imagin g-of-blood-clot-during-microbubble/10.1117/1.JBO.24.12.121902.full
- Mantri Y, Dorobek TR, Tsujimoto J, Penny WF, Garimella PS, Jokerst J V. Monitoring peripheral hemodynamic response to changes in blood pressure via photoacoustic imaging. Photoacoustics [Internet]. 2022; 26: 100345. Available
- https://linkinghub.elsevier.com/retrieve/pii/S2213597922000155
 6. Fadhel MN, Appak Baskoy S, Wang Y, Hysi E, Kolios MC. Use of photoacoustic imaging for monitoring vascular disrupting cancer treatments. J Biophotonics [Internet]. 2023; 16. Available at: https://onlinelibrary.wiley.com/doi/10.1002/jbio.202000209
- Zhu Y, Lu X, Dong X, Yuan J, Fabiilli ML, Wang X. LED-Based Photoacoustic Imaging for Monitoring Angiogenesis in Fibrin Scaffolds. Tissue Eng Part C Methods [Internet]. 2019; 25: 523–31. Available at: https://www.liebertpub.com/doi/10.1089/ten.tec.2019.0151
- Karthikesh MS, Yang X. Photoacoustic image-guided interventions. Exp Biol Med [Internet]. 2020; 245: 330–41. Available at: http://journals.sagepub.com/doi/10.1177/1535370219889323
- Moore C, Chen F, Wang J, Jokerst J V. Listening for the therapeutic window: Advances in drug delivery utilizing photoacoustic imaging. Adv Drug Deliv Rev [Internet]. 2019; 144: 78–89. Available at: https://linkinghub.elsevier.com/retrieve/pii/S0169409X19300845
- Montazar L, Mahani M, Khakbaz F, Divsar F, Yoosefian M. Carbon quantum dots and gold nanostructures on photothermal therapy for cancer treatment. J Photochem Photobiol A Chem [Internet]. 2023; 115443. Available at: https://linkinghub.elsevier.com/retrieve/pii/S1010603023009085
- Shi Y, Zhu D, Wang D, et al. Recent advances of smart AlEgens for photoacoustic imaging and phototherapy. Coord Chem Rev [Internet]. 2022; 471: 214725. Available at: https://linkinghub.elsevier.com/retrieve/pii/S0010854522003204
- Algi MP, Sarıgöl R. Cross-linker engineering of hydrogel enables photothermal therapy and controlled drug release. J Drug Deliv Sci Technol [Internet]. 2023; 89: 104993. Available at: https://linkinghub.elsevier.com/retrieve/pii/S1773224723008456
- Liu P, Mu X, Zhang X-D, Ming D. The Near-Infrared-II Fluorophores and Advanced Microscopy Technologies Development and Application in Bioimaging. Bioconjug Chem [Internet]. 2020; 31: 260–75. Available at: https://pubs.acs.org/doi/10.1021/acs.bioconjchem.9b00610
- Dai H, Shen Q, Shao J, Wang W, Gao F, Dong X. Small Molecular NIR-II Fluorophores for Cancer Phototheranostics. Innov [Internet]. 2021; 2: 100082.

Available

https://linkinghub.elsevier.com/retrieve/pii/S2666675821000072

 Wang J, Wu Y, Liu K, et al. PH/thermal dual-responsive multifunctional drug delivery system for effective photoacoustic imaging-guided tumor chemo/photothermal therapy. APL Bioeng. 2023; 7: 1–11.

at:

- Wang Q, Dai Y, Xu J, et al. All-in-One Phototheranostics: Single Laser Triggers NIR-II Fluorescence/Photoacoustic Imaging Guided Photothermal/Photodynamic/Chemo Combination Therapy. Adv Funct Mater. 2019; 29: 1–12.
- Xu C, Wang Y, Wang E, et al. Effective Eradication of Tumors by Enhancing Photoacoustic-Imaging-Guided Combined Photothermal Therapy and Ultrasonic Therapy. Adv Funct Mater. 2021; 31: 1–10.
- Zhang R, Wang L, Wang X, et al. Acid-Induced *In vivo* Assembly of Gold Nanoparticles for Enhanced Photoacoustic Imaging-Guided Photothermal Therapy of Tumors. Adv Healthc Mater. 2020; 9: 1–10.
- Ye F, Huang W, Li C, et al. Near-Infrared Fluorescence/Photoacoustic Agent with an Intensifying Optical Performance for Imaging-Guided Effective Photothermal Therapy. Adv Ther. 2020; 3: 1–10.
- Zhang Y, Wang Y, Yang X, Yang Q, Li J, Tan W. Polyaniline Nanovesicles for Photoacoustic Imaging-Guided Photothermal-Chemo Synergistic Therapy in the Second Near-Infrared Window. Small. 2020; 16: 1–9.
- Zhou Z, Li B, Shen C, et al. Metallic 1T Phase Enabling MoS2 Nanodots as an Efficient Agent for Photoacoustic Imaging Guided Photothermal Therapy in the Near-Infrared-II Window. Small. 2020; 16: 1–9.
- 22. Chang K, Gao D, Qi Q, Liu Y, Yuan Z. Engineering biocompatible benzodithiophene-based polymer dots with tunable absorptions as high-efficiency theranostic agents for multiscale photoacoustic imaging-guided photothermal therapy. Biomater Sci. 2019; 7: 1486–92.
- Ou C, Zhang Y, Ge W, et al. A three-dimensional BODIPY-iron(III) compound with improved H2O2-response for NIR-II photoacoustic imaging guided chemodynamic/photothermal therapy. Chem Commun. 2020; 56: 6281-4.
- Zheng Z, Chen Q, Dai R, et al. A continuous stimuli-responsive system for NIR-II fluorescence/photoacoustic imaging guided photothermal/gas synergistic therapy. Nanoscale. 2020; 12: 11562–72.
- Sun J, Li Y, Teng Y, Wang S, Guo J, Wang C. NIR-controlled HSP90 inhibitor release from hollow mesoporous nanocarbon for synergistic tumor photothermal therapy guided by photoacoustic imaging. Nanoscale. 2020; 12: 14775–87.
- 26. Sun J, Cai W, Sun Y, Guo C, Zhang R. Facile synthesis of melanin-dye nanoagent for NIR-II fluorescence/photoacoustic imaging-guided photothermal therapy. Int J Nanomedicine [Internet]. 2020; 15: 10199-213. Available at: https://www.dovepress.com/facile-synthesis-of-melanin-dye-nanoagent-forinternet.

nir-ii-fluorescencephoto-peer-reviewed-article-IJN 27. Jun SW, Manivasagan P, Kwon J, et al. Folic acid-conjugated chitocan functionalized graphene oxide for highly efficient photoacoustic

- chitosan-functionalized graphene oxide for highly efficient photoacoustic imaging-guided tumor-targeted photothermal therapy. Int J Biol Macromol [Internet]. 2020; 155: 961–71. Available at: https://doi.org/10.1016/j.ijbiomac.2019.11.055
- Li J, Jiang R, Wang Q, et al. Semiconducting polymer nanotheranostics for NIR-II/Photoacoustic imaging-guided photothermal initiated nitric oxide/photothermal therapy. Biomaterials [Internet]. 2019; 217: 119304. Available at: https://doi.org/10.1016/j.biomaterials.2019.119304
- Wang D, Zhang Z, Lin L, et al. Porphyrin-based covalent organic framework nanoparticles for photoacoustic imaging-guided photodynamic and photothermal combination cancer therapy. Biomaterials [Internet]. 2019; 223: 119459. Available at: https://doi.org/10.1016/j.biomaterials.2019.119459
 Wu J, Williams GR, Zhu Y, et al. Ultrathin chalcogenide nanosheets for
- Wu J, Williams GR, Zhu Y, et al. Ultrathin chalcogenide nanosheets for photoacoustic imaging-guided synergistic photothermal/gas therapy. Biomaterials [Internet]. 2021; 273: 120807. Available at: https://doi.org/10.1016/j.biomaterials.2021.120807
- Rong X, Xia X, Wang R, et al. Near-infrared and highly photostable squaraine-based nanoparticles for photoacoustic imaging guided photothermal therapy. Dye Pigment [Internet]. 2023; 211: 111055. Available at: https://doi.org/10.1016/j.dyepig.2022.111055
- Han Y, Qu B, Li J, et al. A simple POM clusters for *in vivo* NIR-II photoacoustic imaging-guided NIR-II photothermal therapy. J Inorg Biochem [Internet]. 2020; 209: 111121. Available at: https://doi.org/10.1016/j.jinorgbio.2020.111121
- He Y, Pan Y, Zhao X, Fan W, Cai Y, Mou X. NIR-II absorptive dithienopyrrole-thiadiazolobenzotriazole conjugated polymer for photoacoustic imaging-guided glioblastoma multiforme photothermal therapy. Acta Biomater. 2022; 152: 546–61.
- Chen C, Chen Y, Wang X, et al. In situ synthesized nanozyme for photoacoustic-imaging-guided photothermal therapy and tumor hypoxia relief. iScience [Internet]. 2023; 26: 106066. Available at: https://doi.org/10.1016/j.isci.2023.106066
- 35. Zhang R, Wang Z, Xu L, et al. Rational Design of a Multifunctional Molecular Dye with Single Dose and Laser for Efficiency NIR-II Fluorescence/Photoacoustic Imaging Guided Photothermal Therapy. Anal Chem. 2019; 91: 12476–83.
- Xing X, Yang K, Li B, et al. Boron Dipyrromethene-Based Phototheranostics for Near Infrared Fluorescent and Photoacoustic Imaging-Guided Synchronous Photodynamic and Photothermal Therapy of Cancer. J Phys Chem Lett. 2022; 13: 7939–46.

- Guo T, Tang Q, Guo Y, et al. Boron Quantum Dots for Photoacoustic Imaging-Guided Photothermal Therapy. ACS Appl Mater Interfaces. 2021; 13: 306–11.
- Yang K, Long F, Liu W, et al. A-DA'D-A Structured Organic Phototheranostics for NIR-II Fluorescence/Photoacoustic Imaging-Guided Photothermal and Photodynamic Synergistic Therapy. ACS Appl Mater Interfaces. 2022; 14: 18043–52.
- Wang Y, Yang Y, Yang L, et al. Gold Nanostar@Polyaniline Theranostic Agent with High Photothermal Conversion Efficiency for Photoacoustic Imaging-Guided Anticancer Phototherapy at a Low Dosage. ACS Appl Mater Interfaces [Internet]. 2022; 14: 28570–80. Available at: https://pubs.acs.org/doi/10.1021/acsami.2c05679
- 40. Zhang LP, Kang L, Li X, Liu S, Liu T, Zhao Y. Pyrazino[2,3-g]quinoxaline-Based Nanoparticles as Near-Infrared Phototheranostic Agents for Efficient Photoacoustic-Imaging-Guided Photothermal Therapy. ACS Appl Nano Mater. 2021; 4: 2019–29.
- Guo H, Cao Z, Li J, et al. Integrating Bacteria with a Ternary Combination of Photosensitizers for Monochromatic Irradiation-Mediated Photoacoustic Imaging-Guided Synergistic Photothermal Therapy. ACS Nano. 2023; 17: 5059-71.
- 42. Li CQ, Hou XL, Jiang DX, et al. Binary Pt/Te Nanoheterostructures with High Photothermal Conversion Efficiency and Anti-inflammatory Action for Enhanced Photothermal Therapy of 4T1 Breast Tumors Guided by Photoacoustic Imaging. ACS Sustain Chem Eng. 2022; 10: 16598–610.
- Li B, Pang E, Zhao S, et al. Aggregation-Enhanced Photophysical Performance of D-n-A Structured Hemicyanine for NIR-II Fluorescent and Photoacoustic Imaging-Guided Photothermal Therapy. Chem Biomed Imaging. 2023; 1: 541– o
- Chen Y, Wang M, Zheng K, et al. Antimony Nanopolyhedrons with Tunable Localized Surface Plasmon Resonances for Highly Effective Photoacoustic-Imaging-Guided Synergistic Photothermal/Immunotherapy. Adv Mater. 2021; 33: 1–10.
- Yang G, Mu X, Pan X, et al. Ligand engineering of Au44 nanoclusters for NIR-II luminescent and photoacoustic imaging-guided cancer photothermal therapy. Chem Sci. 2023; 14: 4308–18.
- Dai X, Zhao X, Liu Y, et al. Controlled Synthesis and Surface Engineering of Janus Chitosan-Gold Nanoparticles for Photoacoustic Imaging-Guided Synergistic Gene/Photothermal Therapy. Small. 2021; 17: 1–14.
- Wang Z, Sun X, Huang T, Song J, Wang Y. A Sandwich Nanostructure of Gold Nanoparticle Coated Reduced Graphene Oxide for Photoacoustic Imaging-Guided Photothermal Therapy in the Second NIR Window. Front Bioeng Biotechnol. 2020; 8: 1–10.
- Gao F, Miao Y, Ma H, Zhang T, Fan H, Zhao L. Boosting the photothermal performance of vacancy-rich MoSe2-: Xnanoflowers for photoacoustic imaging guided tumor chemo-photothermal therapy. Nanoscale. 2021; 13: 14960-72.
- Guo Q, Wang D, Yang G. Photoacoustic imaging guided photothermal and chemodynamic combined therapy for cancer using 'All in one' W18O49 nanorod agent. J Biomed Nanotechnol. 2019; 15: 2090–9.
- Qiu T, Lan Y, Wei Z, et al. *In vivo* multi-scale photoacoustic imaging guided photothermal therapy of cervical cancer based on customized laser system and targeted nanoparticles. Int J Nanomedicine. 2021; 16: 2879–96.
- Jin X, Xing X, Deng Q, Qing W, Liu Z, Huang Y. Molecular engineering of diketopyrrolopyrrole-conjugated polymer nanoparticles by chalcogenide variation for photoacoustic imaging guided photothermal therapy. J Mater Chem B. 2021; 9: 3153–60.
- Liang G, Han J, Xing D. Precise Tumor Photothermal Therapy Guided and Monitored by Magnetic Resonance/Photoacoustic Imaging using A Safe and pH-Responsive Fe(III) Complex. Adv Healthc Mater. 2021; 10: 1–13.
- Liu G, Zhu J, Guo H, et al. Mo2C-Derived Polyoxometalate for NIR-II Photoacoustic Imaging-Guided Chemodynamic/Photothermal Synergistic Therapy. Angew Chemie - Int Ed. 2019; 58: 18641–6.
- Liu F, Ma F, Chen Q, et al. Synergistic non-bonding interactions based on diketopyrrolo-pyrrole for elevated photoacoustic imaging-guided photothermal therapy. Biomater Sci. 2021; 9: 908–16.
- Zhang L, Yang P, Guo R, Sun J, Xie R, Yang W. Multifunctional mesoporous polydopamine with hydrophobic paclitaxel for photoacoustic imaging-guided chemo-photothermal synergistic therapy. Int J Nanomedicine. 2019; 14: 8647– 63.
- Ou C, Na W, Ge W, et al. Biodegradable Charge-Transfer Complexes for Glutathione Depletion Induced Ferroptosis and NIR-II Photoacoustic Imaging Guided Cancer Photothermal Therapy. Angew Chemie. 2021; 133: 8238–44.
- Li B, Gan Y, Yang K, et al. Acceptor-donor-acceptor structured phototheranostics for near-infrared II fluorescent and photoacoustic imaging-guided photodynamic and photothermal synergistic therapy. Sci China Mater. 2023; 66: 385–94.
- Tian R, Zhu L, Qin Z, Wang G, Wang J, Zhang H. Glypican-3 (GPC3) targeted Fe3O4 core/Au shell nanocomplex for fluorescence/MRI/photoacoustic imaging-guided tumor photothermal therapy. Biomater Sci. 2019; 7: 5258–69.
- Men X, Chen H, Sun C, et al. Thermosensitive Polymer Dot Nanocomposites for Trimodal Computed Tomography/Photoacoustic/Fluorescence Imaging-Guided Synergistic Chemo-Photothermal Therapy. ACS Appl Mater Interfaces. 2020; 12: 51174–84.
- 60. Zhang L, Jiang C, Li B, et al. A core-shell Au@Cu2-xSe heterogeneous metal nanocomposite for photoacoustic and computed tomography

dual-imaging-guided photothermal boosted chemodynamic therapy. J Nanobiotechnology [Internet]. 2021; 19: 1–18. Available at: https://doi.org/10.1186/s12951-021-01159-x

- Li M, Bian X, Chen X, et al. Multifunctional liposome for photoacoustic/ultrasound imaging-guided chemo/photothermal retinoblastoma therapy. Drug Deliv [Internet]. 2022; 29: 519–33. Available at: https://doi.org/10.1080/10717544.2022.2032876
- Lin Y, Zhong W, Wang M, Chen Z, Lu C, Yang H. Multifunctional Carbon Monoxide Prodrug-Loaded Nanoplatforms for Effective Photoacoustic Imaging-Guided Photothermal/Gas Synergistic Therapy. ACS Appl Bio Mater. 2021; 4: 4557–64.
- Chen Y, Xu C, Cheng Y, Cheng Q. Photostability enhancement of silica-coated gold nanostars for photoacoustic imaging guided photothermal therapy. Photoacoustics [Internet]. 2021; 23: 100284. Available at: https://doi.org/10.1016/j.pacs.2021.100284
- Sun L, Jiao X, Liu W, et al. Novel Oxygen-Deficient Zirconia (ZrO2- x) for Fluorescence/Photoacoustic Imaging-Guided Photothermal/Photodynamic Therapy for Cancer. ACS Appl Mater Interfaces. 2019; 11: 41127–39.
- Sun J, Li L, Cai W, Chen A, Zhang R. Multifunctional Hybrid Nanoprobe for Photoacoustic/PET/MR Imaging-Guided Photothermal Therapy of Laryngeal Cancer. ACS Appl Bio Mater. 2021; 4: 5312–23.
- Hu X, Tang Y, Hu Y, et al. Gadolinium-chelated conjugated polymer-based nanotheranostics for photoacoustic/magnetic resonance/NIR-II fluorescence imaging-guided cancer photothermal therapy. Theranostics. 2019; 9: 4168–81.
- Li C, Yang XQ, An J, et al. A near-infrared light-controlled smart nanocarrier with reversible polypeptide-engineered valve for targeted fluorescence-photoacoustic bimodal imaging-guided chemo-photothermal therapy. Theranostics. 2019; 9: 7666–79.
- Zheng T, Zhou T, Feng X, Shen J, Zhang M, Sun Y. Enhanced Plasmon-Induced Resonance Energy Transfer (PIRET)-Mediated Photothermal and Photodynamic Therapy Guided by Photoacoustic and Magnetic Resonance Imaging. ACS Appl Mater Interfaces. 2019; 11: 31615–26.
- Jiang Y, Upputuri PK, Xie C, et al. Broadband Absorbing Semiconducting Polymer Nanoparticles for Photoacoustic Imaging in Second Near-Infrared Window. Nano Lett [Internet]. 2017; 17: 4964–9. Available at: https://pubs.acs.org/doi/10.1021/acs.nanolett.7b02106
- Pulaski BA, Ostrand-Rosenberg S. Mouse 4T1 Breast Tumor Model. Curr Protoc Immunol [Internet]. 2000; 39. Available at: https://currentprotocols.onlinelibrary.wiley.com/doi/10.1002/0471142735.i m2002s39
- Chicheł A, Skowronek J, Kubaszewska M, Kanikowski M. Hyperthermia description of a method and a review of clinical applications. Reports Pract Oncol Radiother [Internet]. 2007; 12: 267–75. Available at: https://linkinghub.elsevier.com/retrieve/pii/S150713671060065X
- Cherukuri P, Glazer ES, Curley SA. Targeted hyperthermia using metal nanoparticles. Adv Drug Deliv Rev [Internet]. 2010; 62: 339–45. Available at: https://linkinghub.elsevier.com/retrieve/pii/S0169409X09003445
- Deng H, Wang L, Tang D, Zhang Y, Zhang L. Review on the laser-induced performance of photothermal materials for ignition application. Energ Mater Front [Internet]. 2021; 2: 201–17. Available at: https://linkinghub.elsevier.com/retrieve/pii/S2666647221000476
- Jain PK, Lee KS, El-Sayed IH, El-Sayed MA. Calculated Absorption and Scattering Properties of Gold Nanoparticles of Different Size, Shape, and Composition: Applications in Biological Imaging and Biomedicine. J Phys Chem B [Internet]. 2006; 110: 7238–48. Available at: https://pubs.acs.org/doi/10.1021/jp057170o